




A Full Range of Products for OR and Reprocessing



LogiCell® Cabinetry in the OR



Case Carts and Par Wall® in the OR



LogiCell® Room Carts in the OR



Basket-Carts and Par Wall® in Reprocessing



Case Carts, Storage and Wrap Carts in Reprocessing

Table of Contents

Case and Reprocessing Carts

2-12



Stainless Steel Furniture and Tables

13-20



LogiCell Stainless Steel Cabinetry

21-27



Work Stations, Sinks & Warming Cabinets

28-37



OR and Reprocessing Storage Systems

38-51



Terms and Conditions / Index

52-53



Case & Reprocessing Carts

Closed Case Carts

LogiQuip Wide Body Case Carts

- Ergonomic by design dual handles.
- New wide shelf, 26", 37".
- Full depth shelf 24".
- Waterfall shelf front.
- Deep for fullsize sets.
- All SS for high temperature cart washing.
- All SS 6" casters, fully cart washable.

Wide Body Closed Carts

MODEL #	OVERALL SIZE D x W x H (inches)	CABINET SIZE D x W x H (inches)
CCWB242630	27.5 x 36 x 30	26.5 x 30.5 x 21.5
CCWB242639	27.5 x 36 x 39	26.5 x 30.5 x 30.25
CCWB242654	27.5 x 36 x 54	26.5 x 30.5 x 45.25
CCWB243739	27.5 x 47 x 39	26.5 x 41.75 x 30.25
CCWB243754	27.5 x 47 x 54	26.5 x 41.75 x 45.25
CCWB244439	27.5 x 54 x 39	26.5 x 48.75 x 30.25
CCWB244454	27.5 x 54 x 54	26.5 x 48.75 x 45.25




Oversized instrument sets side-by-side on shelf

Designed with More Capacity!



Designed for Mobility



CCWB242654 w/ 2426WRS

Case & Reprocessing Carts

Closed Case Carts

Shelf Options for Wide Body Carts

SHELVES	D x W	DESCRIPTION
model #	(inches)	
2426WS	24 x 26	Wire Adjustable for CCWB2426 (30, 39, 54)
2426WSR	24 x 26	Wire Roll-Out for CCWB2426 (30, 39, 54)
2426FS	24 x 26	Solid Adjustable for CCWB2426 (30, 39, 54)
2426FSR	24 x 26	Solid Roll-Out for CCWB2426 (30, 39, 54)
2437WS	24 x 37	Wire Adjustable for CCWB2437 (39, 54)
2437WSR	24 x 37	Wire Roll-Out for CCWB2437 (39, 54)
2437FS	24 x 37	Solid Adjustable for CCWB2437 (39, 54)
2437FSR	24 x 37	Solid Roll-Out for CCWB2437 (39, 54)
2444WS	24 x 44	Wire Adjustable for CCWB2437 (39, 54)
2444WSR	24 x 44	Wire Roll-Out for CCWB2437 (39, 54)
2444FS	24 x 44	Solid Adjustable for CCWB2437 (39, 54)
2444FSR	24 x 44	Solid Roll-Out for CCWB2437 (39, 54)

30" High L & D Cart



CCWB242630

39" High BackTable,
Drapable



CCWB242639



CCWB243739



CCWB244439

54" High Capacity,
Tall Carts



CCWB242654



CCWB243754



CCWB244454

Case & Reprocessing Carts

Closed Case Carts

Industry Standard Size Carts

Standard Size Case Carts

MODEL #	DIMENSIONS OD	CABINET DIMENSIONS	SHELF SIZE	FOOTPRINT	CAPACITY
	D x W x H (inches)	D x W x H (inches)	D x W (inches)	square feet	cubic feet
CCS242430	26.5 x 29.75 x 30	25 x 25 x 21.5	24 x 22.5	4.34	7.78
CCS242439	26.5 x 29.75 x 39	25 x 25 x 30.25	24 x 22.5	4.34	10.94
CCS242454	26.5 x 29.75 x 54	25 x 25 x 45.25	24 x 22.5	4.34	16.37
CCS243639	26.5 x 40.5 x 39	25 x 36 x 30.25	24 x 33.5	7.45	15.76
CCS243654	26.5 x 40.5 x 54	25 x 36 x 45.25	24 x 33.5	7.45	23.57
CCS244839	26.5 x 52.5 x 39	25 x 48 x 30.25	24 x 45.5	9.66	21.01
CCS244854	26.5 x 52.5 x 54	25 x 48 x 45.25	24 x 45.5	9.66	31.42

No Upcharge for Custom Carts

Shelf Options for Standard Carts

SIZE	WIRE ROLL-OUT	WIRE, ADJUST	SOLID SS ROLL-OUT
D x W (inches)	model #	model #	model #
24 x 22.5	2424-WSR	2424-WS	2424-FSR
24 x 33.5	2436-WSR	2436-WS	2436-FSR
24 x 45.5	2448-WSR	2448-WS	2448-FSR



CCS242454



CCS243639

Case & Reprocessing Carts

Closed Case Cart Accessories



Card Holders
8.5" x 11"



Door Seal Safety Tabs



Extendable Drop
Shelf



Tow Hitch

Available Options, Wide Body & Standard

MODEL #	AVAILABLE OPTIONS
CC-CARD	Card holders, 8.5" x 11" stainless
CC-SST	Door safety seal tab, stainless
CC-SDG	Silicon door gasket
CC-TSE	Top extension drop shelf, 16" x 24"
CC-TOW	Tow hitch system, customizable



Silicone Door Gasket

Casters

MODEL #	DESCRIPTION
CCCAST6	Case cart caster, swivel
CCCAST6B	Case cart caster, with brake
CCCAST6DL	Case Cart Caster, directional lock



CCCAST6



CCCAST6B

Standard Casters

- SS, 6" x 1.5".
- Sealed bearings, upper & lower race.
- Grease fitting, bearings and axle.
- High temp wheels, poly.
- Cart washable.
- Standard configuration is all swivel, one directional lock, two braking.
- Optional sizes available upon request.



CCCAST6DL

Case & Reprocessing Carts

Open Case Carts



CCL2436T



CCL2940



CCL2853



OCW2436



Open Case Carts are the best solution for dedicated clean and soiled delivery systems between the OR and reprocessing. Lighter, easier to clean and more cost-effective than closed carts. Same 5-year guarantee, 3-years on casters.

Open Case Carts, 39" Tall

MODEL #	D x W x H (inches)	DESCRIPTION
		cart only
CCL2436T	26 x 36 x 39	Integrated handles on tubular frame, fixed position top/solid shelf, adjustable middle shelf option
CCL2940	29 x 40 x 39	Work top cart with fixed top shelf, adjustable middle shelf optional
CCL2853	28 x 53 x 39	Work top cart with fixed top shelf, adjustable middle shelf optional
CCL2436TFR	29 x 36 x 39	Integrated handles on tubular frame, adjustable top & middle shelves
OCW2436	28 x 40 x 39	All welded, solid SS, fixed position shelf, 6" SS casters, shelf size 24" x 36"



CCL2436TFR
shown w/ optional accessories

Case & Reprocessing Carts

Open Case Carts



CCL293454T
shown w/ optional shelf



CCL293654T
shown w/ optional shelf



CCL266067T
shown w/ optional shelf

Tall Open Case Carts increase the load capacity of instrumentation & supplies.

Open Case Carts, 54" & 67" Tall

MODEL #	OVERALL SIZE	DESCRIPTION
	D x W x H (inches)	Cart Only
CCL293454T	29 x 34 x 54	Large Capacity in verticle sapce. Fully adjustable shelves, push handles built into tubular frame.
CCL293654T	29 x 36 x 54	Extra large capacity for instruments and supplies. Fully adjustable shelves, push handles built into tubular frame.
CCL266067T	26 x 60 x 67	Extra wide cart for very large cases on a 54" wide shelf. Tall design to allow XL custom packs for heart cases.

Open Case Cart Shelves


31" CART	DIMENSIONS	WIRE	SOLID
model #	D x W (inches)	model #	model #
CCL2436T	25.5 x 37	WL2637	SL2637
CCL2940	22 x 36	WL2236	SL2236
CCL2853	21 x 48	WL2148	SL2148
CCL2436TFR	20.75 x 35	WL2135	SL2135
OCW2436	25 x 36	n/a	integrated in

Open Case Cart Shelves

54", 67" CART	DIMENSIONS	WIRE	SOLID
model #	D x W (inches)	model #	model #
CCL293454T	24 x 30	WL2430	SL2430
CCL293654T	24 x 32	WL2432	SL2432
CCL266067T	24 x 56	WL2456	SL2456

Case & Reprocessing Carts

Open Case Cart Accessories



Backscreen of 3" x 3" wire mesh welded to frame.



Card Holders



Cart washable baskets and totes



Different caster heights and configurations



Nylon or disposable plastic cart covers (clear or opaque).

Case & Reprocessing Carts

The Over the Road Case Cart

Total Stainless Cart Transport sterile supplies and instrumentation offsite safely and cleanly.

- 30"D x 43.75"W x 57"H OD (8.28 sq. ft.)
- 28"D x 39.25"W x 45.75"H cabinet (24.94 cu. ft.)
- Shelf size: 26"D x 39"W.
- Capacity for 6 full size containers.
- 12, 18, 24 container capacity available.
- Unique container retention system.
- Top storage space: 5"H.
- Bottom storage space: 14.5"H.
- All welded stainless construction.
- 6" x 2" heavy-duty cart washable casters.
- Two directional locking, two swivel with brake.
- Full wrap-around bumper system on top and bottom of cart.
- Overlapping self-sealing door system, with 95° cart seal.
- Double locking key doors.
- Breakaway security tab.

Model #OTRCC3657SS

1/2 size cart Model #OTRCC2657SS

30"D x 31.5"W x 57"H OD



Easy slide indicator eliminates
guesswork and errors



OTRCC3657SS



OTRCC3657SS



Flag System



Tow Hitch



OTRCC2657SS

Case & Reprocessing Carts

Wash Carts



UW101-4HP

Wash Carts

- All adjustable for varied size containers
- All welded, type 304 stainless
- Hi temp 5" phenolic casters, SS with grease fitting
- High capacity
- Ergonomic design
- Open design for full wash of containers, pans, lids, and filters

MODEL #	DIMENSIONS	MAX CAPACITY	DESCRIPTION
	D x W x H (inches)		
UW101-4HP	25 x 60.5 x 69	16	Soft coat rack
UW101-4WST-61	21 x 61.5 x 55	16	All welded construction. Capacity design for containers, pans and misc reusables
UW101-4WST-80	21 x 80.5 x 55	20	All welded construction. Capacity design for containers, pans and misc reusables



UW101-4WST61



UW101-4WST80

Case & Reprocessing Carts

Wrap/Pick Carts



WC60CH

Sterile Wrap Rack

- Chrome plated steel.
- 60" models hold wraps up to 54" x 72".
- Four bars standard. Rubber "O" rings attached to bars to secure wraps from sliding off rack.
- 5" swivel non-marring casters, 3" non-marring bumpers.
- Bottom chrome wire shelf w/ plastic shelf cover.

Model #WC48CH: 24"D x 48"W x 69"H.

Model #WC60CH: 24"D x 60"W x 69"H.

Assembly Cart

- For use beside every reprocessing worktable.
- 24"D x 48"W x 80"H.
- Five chrome wire shelves.
- Solid plastic bottom shelf cover.
- 5" swivel non-marring casters, two with brakes.

Model #MB455CH-CL



MB455CH-CL

Express/Pick Carts

- Epoxy-Coated Steel Fabrication.
- Ergonomically designed push handle at 39-1/2".
- Two fixed, two swivel 5" casters.
- Standard color: Gray.



Express/Pick Carts Two Baskets

- Overall Dimensions: 27-1/2"D x 20-1/2"W x 39-1/2"H.
- Top Basket Dimensions: 17-1/2"D x 16-1/4"W x 8"H.
- Bottom Basket Dimensions: 17-1/2"D x 16-1/4"W x 4"H.
- Weight: 39 lbs.

Model #5341



Express/Pick Carts One Basket

- Overall Dimensions: 22-1/2"D x 20"W x 39-1/2"H.
- Basket Dimensions: 17-1/2"D x 16-1/4"W (14"W at front) x 8"H.
- Weight: 30 lbs.

Model #5141

Case & Reprocessing Carts

Ortho Cart System



6200-III Container Cart

Ortho Cart System Features and Benefits

- A clean cart system designed to store and transport various types of ortho implants and/or related instrument sets within surgical services.
- Carts are maintained in reprocessing and sent to the OR with a companion case cart. The case cart holds basic sets and picked supplies.
- Eliminates/reduces the need for giant ortho case carts or multiple smaller case carts.
- Integrates well with loaner instrument programs.



6200-III Cart with Implants

Ortho Implant and Container Cart Specifications

- All-welded 1-1/2" square tubular epoxy-coated steel frame.
- All shelves waterfall-edge and baskets on full-extension slides. Stoppers available.
- Corner bumpers, 4 x 5" swivel casters, two with brakes.

MODEL #	DESCRIPTION	D x W x H (O.D.) (inches)
3000-III	Single short implant cart. 1 - 3", 2 - 5", 1 - 8" baskets. Optional SS top available	27.5 x 25 x 43
PS300-III	Single short container cart. 3 chrome wire shelves & SS top current build has 3 shelves & a solid top	27.5 x 25 x 43
3200-III	Double short implant cart. 2 - 3", 4 - 5", 2 - 8" baskets. Optional SS top available	27.5 x 44 x 43
PS320-III	Double short container cart. 8 chrome wire shelves. Optional SS top available	27.5 x 44 x 43
5000-III	Single medium implant cart. 1 - 3", 4 - 5", 1 - 8" baskets	27.5 x 25 x 58.5
PS500-III	Single medium container cart. 6 chrome wire shelves	27.5 x 25 x 58.5
5200-III	Double medium implant cart. 2 - 3", 8 - 5", 2 - 8" baskets	27.5 x 44 x 58.5
PS520-III	Double medium container cart. 12 chrome wire shelves	27.5 x 44 x 58.5
6000-III	Single tall implant cart. 1 - 3", 6 - 5" baskets	27.5 x 25 x 74.5
PS600-III	Single tall container cart. 7 chrome wire shelves	27.5 x 25 x 74.5
6200-III	Double tall implant/container cart. 7 chrome wire shelves and 7 - 5" baskets	27.5 x 44 x 74.5
PS620-III	Double tall container cart. 14 each chrome wire shelves	27.5 x 44 x 74.5
6300-III	Triple tall implant/container cart. 7 chrome wire shelves and 14 - 5" baskets	27.5 x 65 x 74.5
PS630-III	Triple tall container cart. 21 each chrome wire shelves	27.5 x 65 x 74.5

Not automatic cart washable, hand wash only.

Stainless Steel Furniture

Back - Instrument Tables

MODEL #	D x W x H	DESCRIPTION
	(inches)	
BT1620	16 x 20 x 34	SS back/instrument table, bottom shelf
BT1630	16 x 30 x 34	SS back/instrument table, bottom shelf
BT2436	24 x 36 x 34	SS back/instrument table, bottom shelf
BT2448	24 x 48 x 34	SS back/instrument table, bottom shelf
BT2460	24 x 60 x 34	SS back/instrument table, bottom shelf
BT2472	24 x 72 x 34	SS back/instrument table, bottom shelf
BT1620H	16 x 20 x 34	SS back/instrument table, h-brace
BT1630H	16 x 30 x 34	SS back/instrument table, h-brace
BT2436H	24 x 36 x 34	SS back/instrument table, h-brace
BT2448H	24 x 48 x 34	SS back/instrument table, h-brace
BT2460H	24 x 60 x 34	SS back/instrument table, h-brace
BT2472H	24 x 72 x 34	SS back/instrument table, h-brace

*Additional sizes available for specific requirements.
Please contact us for further assistance.*

Standard Features:

- Fully-welded construction.
- Reinforced tops.
- Durable 1 1/4" 16 gauge stainless steel legs.
- 16 gauge type 304 stainless steel top.
- 18 gauge type 304 stainless steel under shelf.
- 3" swivel casters.

Options:

- Swivel casters.
- Locking casters.
- Drawers mount under worktop.
- Stainless steel guard rails.
- Nesting sets. (see page 15)
- Multiple under shelves.



BT1620



BT1630H



BT2448

Stainless Steel Furniture

Utility Tables - Prep Stands



UT1636D Single Drawer



UT1636D2 Double Drawer



UT1636D3 Triple Drawer

MODEL #	D x W x H	DESCRIPTION
	(inches)	
UT1636D	16 x 20 x 34	SS back/instrument table, 3-sided top rail, 1 drawer, 2 shelves
UT1636D2	16 x 20 x 34	SS back/instrument table, 3-sided top rail, 2 drawers, 2 shelves
UT1636D3	16 x 20 x 34	SS back/instrument table, 3-sided top rail, 3 drawer, 2 shelves
BT1620-GR	16 x 20 x 34	Utility table, 4 sided guardrail on top and bottom

Additional sizes available for specific requirements.

Prep Stands / Utility Tables

- Available in 3 models: single, double and triple drawers.
- Size: 16"D x 20"W x 34"H.
- Caster, 2", 2 with brake.
- Drawers, 4"H.

Stainless Steel Furniture

Work Tables + Nesting Tables

Heavy Duty Work Tables

MODEL #	D x W x H	DESCRIPTION
	(inches)	
MTF 101	30 x 60 x 34	SS work table, H-brace
MTF 102	36 x 60 x 34	SS work table, H-brace
MTF 103	44 x 60 x 34	SS work table, H-brace
MTF 104	30 x 72 x 34	SS work table, H-brace
MTF 105	36 x 72 x 34	SS work table, H-brace
MTF 106	44 x 72 x 34	SS work table, H-brace
MTF 107	30 x 96 x 34	SS work table, H-brace
MTF 108	36 x 96 x 34	SS work table, H-brace
MTF 109	44 x 96 x 34	SS work table, H-brace
MTS 101	30 x 60 x 34	SS work table, bottom shelf
MTS 102	36 x 60 x 34	SS work table, bottom shelf
MTS 103	44 x 60 x 34	SS work table, bottom shelf
MTS 104	30 x 72 x 34	SS work table, bottom shelf
MTS 105	36 x 72 x 34	SS work table, bottom shelf
MTS 106	44 x 72 x 34	SS work table, bottom shelf
MTS 107	30 x 96 x 34	SS work table, bottom shelf
MTS 108	36 x 96 x 34	SS work table, bottom shelf
MTS 109	44 x 96 x 34	SS work table, bottom shelf



MTF101



MTS104

Nesting Back - Instrument Tables

MODEL #	D x W x H	DESCRIPTION
	(inches)	
NT1428	14 x 28 x 32	SS nesting table, U-brace
NT1632	16 x 32 x 34	SS nesting table, U-brace
NT1836	18 x 36 x 36	SS nesting table, U-brace
NT2040	20 x 40 x 38	SS nesting table, U-brace
NT2244	22 x 44 x 40	SS nesting table, U-brace
NT2448	24 x 48 x 42	SS nesting table, U-brace
NT2448-6	24 x 48 x 42	SS nesting table, U-brace, set of 6





Stainless Steel Furniture

Ring Stands + Kick Buckets



DBS



DBSS

Ring Stands

- Available in single and double basin designs.
- Fully welded stainless steel construction.
- Available with lower shelf or H-brace.
- 2" casters swivel.
- Accessible from any side with open frame design.



SBSS



SBS

MODEL #	SIZE	D x W x H	DESCRIPTION
SBS	Single	14 x 14 x 34	SS basin stand, H-brace
SBSS	Single	14 x 14 x 34	SS basin stand, bottom shelf
DBS	Double	14 x 28 x 34	SS basin stand, H-brace
DBSS	Double	14 x 28 x 34	SS basin stand, bottom shelf
SBS-NB	Single	14 x 14 x 34	SS basin stand, no basin, H-brace
SBSS-NB	Single	14 x 14 x 34	SS basin stand, no basin, bottom shelf
DBS-NB	Double	14 x 28 x 34	SS basin stand, no basin, H-brace
DBSS-NB	Double	14 x 28 x 34	SS basin stand, no basin, bottom shelf
SSBASIN	7 Quart		SS basin
SSBASIN-8.5	8.5 Quart		SS basin

Kick Buckets

- Stainless steel frame.
- 2" swivel casters.
- Perimeter bumper to protect walls.
- Tip resistant.



KBSB



KB

Model# KB

- Overall size: 14"H x 13 1/2" diameter.
- Removable 12 quart stainless steel basin.

Model# KBSB

- Overall size: 8"H x 16 1/4" diameter.
- Removable 8 1/2 quart stainless steel basin.

Stainless Steel Furniture

Foot Operated

Double Column Support



MD1621F

Single Column Support



M1621F

Foot Operated, Friction Lock


MODEL #	TRAY SIZE (inches)	COLUMNS # of
M1319F	13.625 x 19.125	Single
MD1319F	13.625 x 19.125	Dual
M1621F	16.25 x 21.25	Single
MD1621F	16.25 x 21.25	Dual
MD2025F	20 x 25	Dual

Hand Operated

MODEL #	TRAY SIZE (inches)	COLUMNS # of	LOCK DESCRIPTION
M1319TK	13.625 x 19.125	Single	Hand operated, knob operated
MD1319H	12.625 x 19.125	Dual	Hand operated, friction lock
M1621TK	16.25 x 21.25	Single	Hand operated, knob operated
MD1621H	16.25 x 21.26	Dual	Hand operated, friction lock
MD2025H	20 x 25	Dual	Hand operated, friction lock
M1319T	12.625 x 19.125	Single	Thumb operated, friction lock
M1621T	16.25 x 21.25	Single	Thumb operated, friction lock

Manual

HD Dual Column Support



MD1621H

Single Column Support



M1621TK

Single Column Support



M1621T

Stainless Steel Furniture

Step/Foot Stools + Seating

Step / Foot Stools

Model# 121608C: Chrome Step Stool with Black Rubber Mat

- 12"W x 16"L x 8-1/4"H.
- 1" 16-gauge tubular steel legs with rubber floor tips.
- Slip-resistant surface and tip-resistant design.

Model# 121608C-HR: Same as above but with 34"H steel handrail.

Model# 121608S: Same as 121608C but stainless steel.

Model# 121608S-HR: Same as above but with 34"H steel handrail.

Model# 14146-STAK: Stackable Stainless Step Stool

- 14"W x 19"L x 6"H- Largest platform available.
- 1" 16-gauge angle steel legs with rubber floor tips.
- Can stack safely for increased height and interlock to create larger platform.

Model# 14146-STAKD: Stackable Step Stool Dolly

- Can transport up to 6 stools. 4 x 3" swivel casters.



11149C
(12608S looks the same except it is stainless steel.)



12608S-HR



14146-STAK
Stackable Footstool

21340B



21521B



Seating

Model #21340B: Surgeon's Stool

- Comfortable Ergonomic Design.
- Hands free, foot operated pneumatic height control.
- 18"-23" Height Adjustment.
- 18" Contour Poly-Foam Cushion.
- 23" Powder coat black aluminum base.
- High impact black plastic backrest shroud.
- Seamless upholstery standard.
- Available with locking casters.
- Available in 50 standard colors.

Model# 21340:

Same as above without back.

Model# 21521B:

Higher 24"-29.5" height adjustment with round foot ring.

Model# 21520B:

Higher 24"-29.5" height adjustment with round foot ring and glides (no casters).

Stainless Steel Furniture

IV Stands


4-Leg Base Units

Hand

MODEL #	BASE	HOOKS # of	DESCRIPTION
LQIVH-2002	4-Leg	2-Hook	IV Stand, Hand Operated
LQIVH-2004	4-Leg	4-Hook	IV Stand, Hand Operated
LQIVH-2006	4-Leg	6-Hook	IV Stand, Hand Operated

Foot

MODEL #	BASE	HOOKS # of	DESCRIPTION
LQIVF-4002	4-Leg	2-Hook	IV Stand, Foot Operated
LQIVF-4004	4-Leg	4-Hook	IV Stand, Foot Operated
LQIVF-4006	4-Leg	6-Hook	IV Stand, Foot Operated



LQIVF-4006

LQIVH-2002


5-Leg Base Units

Hand

MODEL #	BASE	HOOKS # of	DESCRIPTION
LQIVH-1002	5-Leg	2-Hook	IV Stand, Hand Operated
LQIVH-1004	5-Leg	4-Hook	IV Stand, Hand Operated
LQIVH-1006	5-Leg	6-Hook	IV Stand, Hand Operated

Foot

MODEL #	BASE	HOOKS # of	DESCRIPTION
LQIVF-3002	5-Leg	2-Hook	IV Stand, Foot Operated
LQIVF-3004	5-Leg	4-Hook	IV Stand, Foot Operated
LQIVF-3006	5-Leg	6-Hook	IV Stand, Foot Operated



LQIVH-1004

LQIVF-3006

Stainless Steel Furniture

Linen Hampers

1000 Series



Model# LQH-1000

- Overall size: 35"H.
- Accommodates 18" diameter bags.



Model# LQH-1001

- Overall size: 35"H.
- Accommodates 25" diameter bags.

2000 Series



Model# LQH-2000

- Overall size: 35"H.
- Accommodates 18" diameter bag.



Model# LQH-2001

- Overall size: 35"H.
- Accommodates 25" diameter bags.

Stainless Hamper Available Options

- FOL = Foot operated stainless steel lid.
- BG = Rotating bumper guards over casters.
- Capitalize foot and rotating.
- Fully-welded, heavy gauge all stainless steel construction.
- 3" double ball bearing, swivel casters.
- Full, heavy gauge stainless steel platform to support bags & ease of removal.



Model# #HS50

Chrome Hamper Step-On Covered Single Hamper System

Complete System for Identification, Containment and Collection.

- Chrome plated steel, 1" tubing.
- 20.5" x 19" x 38".
- Step-on opening feature.
- 2" casters provide complete mobility.
- Bottom rack to support heavy loads.
- Storage rack for extra bags.

Model #HS50



Model# #CA20

Basic Spring Hamper Stand

- Chrome.
- 12" x 19" x 38".
- 3" casters.
- Plastic lid w/ chrome handle.

Model #CA20 (no lid)

Model # CA21 (with lid)

LogiCell® Stainless Steel Cabinetry

Stainless Steel Cabinetry



40000 Series Stainless Steel Cabinets
with LogiCell cell system for OR, Cath Lab, and IR

A New Concept In OR Storage Design

Stainless Cabinets for OR, Reprocessing, New Hybrid OR's.

- Floor mounted or **mobile**.
- Built in or pass through.
- Full range of locking & entry identification options.
- Standard flat shelf.
- Modular, high density LogiCell system built into cabinets.
- Scope cabinets.
- Catheter cabinets.



40000 Series Mobile
Custom Designs

LogiCell® Stainless Steel Cabinetry

Stainless Steel Moveable Cabinetry with Built-in Cell System

The "No Installation Cost"

Moveable OR, IR, Cath Lab Cabinet System

MODEL #	D x W x H (inches)	DESCRIPTION
22381SSLG	18 x 27 x 82.5	Single, 23 cells, left glass door
22381SSRG	18 x 27 x 82.5	Single, 23 cells, right glass door
22581SSHGLC	18 x 36 x 82.5	Single, 20 cells left, (5) 8" shelves right, 5 dividers, glass doors
22581SSHGRC	18 x 36 x 82.5	Single, 20 cells right, (5) 8" shelves left, 5 dividers, glass doors
22681SSHGLC	18 x 42 x 82.5	Single, 20 cells left, (5) 14" shelves right, 7 dividers, glass doors
22681SSHGRC	18 x 42 x 82.5	Single, 20 cells right, (5) 14" shelves left, 7 dividers, glass doors
24181SSLG	26 x 18 x 82.5	Single, 23 cells, left glass door
24181SSRG	26 x 18 x 82.5	Single, 23 cells, right glass door
24581SSHG	26 x 36 x 82.5	Single, 23 cells, full width, single, glass doors
24582SSHG	26 x 36 x 82.5	Double, 46 cells, glass doors
24883SSHG	26 x 54 x 82.5	Triple, 69 cells, glass doors

Features:

- Roll in stock ready for use.
- Full stainless construction.
- Mobile system, full bumper or internal dolly base.
- Sloped tops.
- Hi-density, modular LogiCell® cell system.

Optional Accessories:

- SS solid shelves, or epoxy coated.
- Wire basket 3", 5", 8" deep.
- Plastic tubs 3", 5", 8" deep.
- Full divider systems.
- Catheter racks.
- Slanted shelves.
- Suture system.



24883SSHG
Full Wraparound Bumper System



24582SSHG integrated dolly base
w/ optional cath slides and shelves

LogiCell® Stainless Steel Cabinetry

Stainless Steel Cabinets



Sloped Top Option



Flat Pan Base

SS Tall Cabinets, Cell System Cabinets

MODEL #	D x W x H	DESCRIPTION
	(inches)	
42381LG	18 x 27 x 84	Cabinet, single column, cells, left glass door
42381RG	18 x 27 x 84	Cabinet, single column, cells, right glass door
42581HG-LC	18 x 36 x 84	Cabinet, cells left, 5 shelves right, 5 shelf dividers, glass doors
42581HG-RC	18 x 36 x 84	Cabinet, cells right, 5 shelves left, 5 shelf dividers, glass doors
42681HG-LC	18 x 42 x 84	Cabinet, cells left, 5 shelves right, 7 shelf dividers, glass doors
42681HG-RC	18 x 42 x 84	Cabinet, cells right, 5 shelves left, 7 shelf dividers, glass doors
44181LG	26 x 19.5 x 84	Cabinet, single column, cells, left glass door
44181RG	26 x 19.5 x 84	Cabinet, single column, cells, right glass door
44582HG	26 x 36 x 84	Cabinet, double column, cells, glass doors

SS Tall Cabinets, Catheter Cabinets

MODEL #	D x W x H	DESCRIPTION
	(inches)	
42580HGCAT1	18 x 36 x 84	Cath cabinet, (7) - 6 hook cath slides, 5 dividers, glass doors
42680HGCAT1	18 x 42 x 84	Cath cabinet, (9) - 6 hook cath slides, 7 dividers, glass doors
44580HGCAT1	26 x 36 x 84	Cath cabinet, (7) - 9 hook cath slides, 5 dividers, glass doors
42580HGCAT2	18 x 36 x 84	Cath cabinet, (7) - 6 hook cath slides
42680HGCAT2	18 x 42 x 84	Cath cabinet, (9) - 6 hook cath slides
44580HGCAT2	26 x 36 x 84	Cath cabinet, (7) - 9 hook cath slides
42580HGCAT3	18 x 36 x 84	Cath cabinet, (14) - 6 hook cath slides, 2 shelves, 4 dividers, glass doors
42680HGCAT3	18 x 42 x 84	Cath cabinet, (18) - 6 hook cath slides, 2 shelves, 4 dividers, glass doors
44580HGCAT3	26 x 36 x 84	Cath cabinet, (14) - 9 hook cath slides, 2 shelves, 4 dividers, glass doors
42580HGCATSH	18 x 36 x 84	Bxd balloon/guide wire cabinet, 5 adj shelves, 25 dividers, glass doors
42680HGCATSH	18 x 42 x 84	Bxd balloon/guide wire cabinet, 5 adj shelves, 35 dividers, glass doors
42780HGCATSH	18 x 48 x 84	Bxd balloon/guide wire cabinet, 5 adj shelves, 45 dividers, glass doors
44580HGCATSH	26 x 36 x 84	Bxd balloon/guide wire cabinet, 5 adj shelves, 25 dividers, glass doors
44680HGCATSH	26 x 42 x 84	Bxd balloon/guide wire cabinet, 5 adj shelves, 35 dividers, glass doors
42580HGCATB	18 x 36 x 84	Bxd cath cabinet, roll-out cath box, top shelf, 5 dividers, glass doors
44580HGCATB	26 x 36 x 84	Bxd cath cabinet, roll-out cath box, top shelf, 5 dividers, glass doors

LogiCell® Stainless Steel Cabinetry

Stainless Steel Cabinets



Tall Cabinets, Flat Pan Base, 5 Adjustable Shelves

64" Height

MODEL #	D x W x H	DESCRIPTION
	(inches)	
42460HG	18 x 30 x 64	Glass doors
42560HG	18 x 36 x 64	Glass doors
42760HG	18 x 48 x 64	Glass doors
43460HG	24 x 30 x 64	Glass doors
43560HG	24 x 36 x 64	Glass doors
43760HG	24 x 48 x 64	Glass doors

Tall Cabinets, Flat Pan Base, 5 Adjustable Shelves

72" Height

MODEL #	D x W x H	DESCRIPTION
	(inches)	
42270LG	18 x 24 x 72	Left glass door
42270RG	18 x 24 x 72	Right glass door
42570HG	18 x 36 x 72	Glass doors
42770HG	18 x 48 x 72	Glass doors
43270LG	24 x 24 x 72	Left glass door
43270RG	24 x 24 x 72	Right glass door
43570HG	24 x 36 x 72	Glass doors
43770HG	24 x 48 x 72	Glass doors

Tall Cabinets, Flat Pan Base, 5 Adjustable Shelves

84" Height

MODEL #	D x W x H	DESCRIPTION
	(inches)	
42280LG	18 x 24 x 84	Left glass door
42280RG	18 x 24 x 84	Right glass door
42380LG	18 x 27 x 84	Left glass door
42380RG	18 x 27 x 84	Right glass door
42580HG	18 x 36 x 84	Glass doors
42680HG	18 x 42 x 84	Glass doors
42780HG	18 x 48 x 84	Glass doors
43280LG	24 x 24 x 84	Left glass door
43280RG	24 x 24 x 84	Right glass door
43580HG	24 x 36 x 84	Glass doors
43780HG	24 x 48 x 84	Glass doors
44180LG	26 x 18 x 84	Left glass door
44180RG	26 x 18 x 84	Right glass door
44280LG	26 x 24 x 84	Left glass door
44280RG	26 x 24 x 84	Right glass door
44580HG	26 x 36 x 84	Glass doors
44680HG	26 x 42 x 84	Glass doors
44780HG	26 x 48 x 84	Glass doors

LogiCell® Stainless Steel Cabinetry

Stainless Steel Cabinets



SS Floor Cabinets with Cell System

MODEL #	D x W x H	DESCRIPTION
	(inches)	
44131LS	26 x 19.5 x 33	Single column, cells, left solid door
44131RS	26 x 19.5 x 33	Single column, cells, right solid door
44532S	26 x 36 x 33	Double column, cells, solid doors

SS Floor Cabinets, Sink Cabinets

MODEL #	D x W x H	DESCRIPTION
	(inches)	
44230S-SNK	26 x 24 x 33	No shelf, solid doors, basin not included
44430S-SNK	26 x 30 x 33	No shelf, solid doors, basin not included
44530S-SNK	26 x 36 x 33	No shelf, solid doors, basin not included

SS Floor Cabinets, 1 Adjustable Shelf

MODEL #	D x W x H	DESCRIPTION
	(inches)	
44130LS	26 x 18 x 33	1 adjustable shelf, left solid door
44130RS	26 x 18 x 33	1 adjustable shelf, right solid door
44230LS	26 x 24 x 33	1 adjustable shelf, left solid door
44230RS	26 x 24 x 33	1 adjustable shelf, right solid door
44430S	26 x 30 x 33	1 adjustable shelf, solid doors
44530S	26 x 36 x 33	1 adjustable shelf, solid doors
44630S	26 x 42 x 33	1 adjustable shelf, solid doors
44730S	26 x 48 x 33	1 adjustable shelf, solid doors

SS Wall Cabinets, 1 Adjustable Shelf

MODEL #	D x W x H	DESCRIPTION
	(inches)	
41120LG	13 x 18 x 30	Left glass door
41120RG	13 x 18 x 30	Right glass door
41220LG	13 x 24 x 30	Left glass door
41220RG	13 x 24 x 30	Right glass door
41420HG	13 x 30 x 30	Glass doors
41520HG	13 x 36 x 30	Glass doors
41620HG	13 x 42 x 30	Glass doors
41720HG	13 x 48 x 30	Glass doors

SS Floor Cabinets, Drawer Cabinets

MODEL #	D x W x H	DESCRIPTION
	(inches)	
44130LS-1D	26 x 18 x 33	1 drawer, 1 adjustable shelf, left solid door
44130RS-1D	26 x 18 x 33	1 drawer, 1 adjustable shelf, right solid door
44230LS-1D	26 x 24 x 33	1 drawer, 1 adjustable shelf, left solid door
44230RS-1D	26 x 24 x 33	1 drawer, 1 adjustable shelf, right solid door
44130-4D	26 x 18 x 33	4 drawers
44230-4D	26 x 24 x 33	4 drawers



LogiCell® Stainless Steel Cabinetry

Nurse Documentation

4000NDWL
Wall Cabinet

4000NDTL
Tall Cabinet

4000NDBL
Base Cabinet

- All stainless construction.
- Sold as 1 unit or individual parts.
- Customizable.
- Vented for computer air circulation.
- Open back for power & IT outlets.



4000NDSLKBFT Complete Unit, Left Hand Orientation

SS Nurse Doc Station Kick Base, Flat Top, Left Hand

MODEL #	D x W x H	DESCRIPTION
	(inches)	
4000NDSL	28 x 70 x 76	Complete unit
4000NDTL	28 x 27 x 76	Tall cabinet, for 4000NDSL
4000NDBL	28 x 43 x 36	Base cabinet, SS worksurface, 4" back splash, for 4000NDSL(-ST)
4000NDWL	17 x 43 x 27	Wall cabinet, for 4000NDSL

SS Nurse Doc Station Kick Base, Slope Top, Right Hand

MODEL #	D x W x H	DESCRIPTION
	(inches)	
4000NDSR	28 x 70 x 76	Complete unit
4000NDTR	28 x 27 x 76	Tall cabinet, for 4000NDSR
4000NDBR	28 x 43 x 36	Base cabinet, SS worksurface, 4" back splash, for 4000NDSR(-ST)
4000NDWR	17 x 43 x 27	Wall cabinet, for 4000NDSR

SS Nurse Doc Station Kick Base, Slope Top, Left Hand

MODEL #	D x W x H	DESCRIPTION
	(inches)	
4000NDSL-ST	28 x 70 x 82	Complete unit
4000NDTL-ST	28 x 27 x 82	Tall cabinet, for 4000NDSL-ST
4000NDBL	28 x 43 x 36	Base cabinet, SS worksurface, 4" back splash, for 4000NDSL(-ST)
4000NDWL-ST	17 x 43 x 33	Wall cabinet, for 4000NDSL-ST

SS Nurse Doc Station Kick Base, Flat Top, Right Hand

MODEL #	D x W x H	DESCRIPTION
	(inches)	
4000NDSR-ST	28 x 70 x 82	Complete unit
4000NDTR-ST	28 x 27 x 82	Tall cabinet, for 4000NDSR-ST
4000NDBR	28 x 43 x 36	Base cabinet, SS worksurface, 4" back splash, for 4000NDSR(-ST)
4000NDWR-ST	17 x 43 x 33	Wall cabinet, for 4000NDSR-ST

LogiCell® Stainless Steel Cabinetry

Scope Cabinets



11636SC9



32636SCHG9



38636SCHG9



LogiCell® Scope Cabinets

- Full body scope racks to protect expensive equipment bodies.
- Fully padded lower cabinet back wall, antimicrobial.
- Drip tray, removable and cleanable.
- Fully vented cabinets.
- Extra capacity over competitors 9 scopes in 36" vs. 8 scopes

Endoscopes Cabinets

SCOPES	D x W x H	TAMBOUR DOOR	GLASS DOOR	STAINLESS GLASS DOOR
Quantity	(inches)	Model #	Model # Right/Left	Model # Right/Left2
6	18 x 27 x 92	11624SC-6	32624SCRG-6 / 32624SCLG-6	42390SCRG-6 / 42390CLG-6
12, 2 levels	18 x 27 x 92	11624SC-12	32624SCRG-12 / 32624SCLG-12	42390SCRG-12 / 42390SCLG-12
9	18 x 36 x 92	11636SC-9	32636SCHG-9	42590SCHG-9
18, 2 levels	18 x 36 x 92	11636SC-18	32636SCHG-18	42590SCHG-18
11	18 x 42 x 92	11642SC-11	32642SCHG-11	42690SCHG-11
22, 2 levels	18 x 42 x 92	11642SC-22	32642SCHG-22	42690SCHG-22

Bronchoscopes Cabinets

SCOPES	D x W x H	TAMBOUR DOOR	GLASS DOOR	STAINLESS GLASS DOOR
Quantity	(inches)	Model #	Model # Right/Left	Model # Right/Left2
5	18 x 27 x 60	11624SC-5	32624SCRG-5 / 32624SCLG-5	42390SCRG-5 / 42390SCLG-5
10, 2 levels	18 x 27 x 60	11624SC-10	32624SCRG-10 / 32624SCLG-10	42390SCRG-10 / 42390SCLG-10
7	18 x 36 x 60	11636SC-7	32636SCHG-7	42590SCHG-7
14, 2 levels	18 x 36 x 60	11636SC-14	32636SCHG-14	42590SCHG-14
9	18 x 42 x 60	11642SC-9	32642SCHG-9	42690SCHG-9
18, 2 levels	18 x 42 x 60	11642SC-18	32642SCHG-18	42690SCHG-18

TEE Probe Cabinets

SCOPES	D x W x H	TAMBOUR DOOR	GLASS DOOR	STAINLESS GLASS DOOR
	(inches)	Model #	Model # Right/Left	Model # Right/Left2
3	18 x 27 x 92	11624SC-3T	32624SCRG-3T / 32624SCLG-3T	42390SCRG-3T / 42390SCLG-3T
5	18 x 36 x 92	11636SC-5T	32636SCHG-5T	42590SCHG-5T
7	18 x 42 x 92	11642SC-7T	32642SCHG-7T	42690SCHG-7T



Double Scope Rack

Sinks, Warming Cabinets & Workstations

SS Series Surgical Scrub Sinks Features and Options



Optional Recommended Products

SINK MODEL	DESCRIPTION	IN-WALL CHAIR CARRIER	SHELF OVER SINK	LIGHTED SHELF OVER SINK
Model #		Model #	Model #	Model #
SS32	Single Sink	WCSS32	SH0001	LSH0001
SS64	Dual Sink	WCSS64	SH0002	LSH0002
SS96	Triple Sink	WCSS96	SH0003	LSH0003

How To Order SS Series Sinks

INFRARED OPERATION	MANUAL BYPASS TO ADD KNEE OPERATED WATER AND SOAP WITH IR OPERATION	KNEE OPERATED SOAP WITH IR OPERATION	FOOT OPERATED SOAP (ALL SINKS)	DIGITAL SCRUB TIMER	EYE WASH	220-240V	PEDESTAL BASE
Opt Ext	Opt Ext	Opt Ext	Opt Ext	Opt Ext	Opt Ext	Opt Ext	Opt Ext
IR	MB	KSP	FSP	TM	EYE	220	PB

For Example: Model # SS64 with the *Infrared Operated with Knee Operated Water and Soap plus Digital Scrub Timers*



Order Model # SS64-IR-MB-TM

Tables, Sinks, and Warming Cabinets

Surgical Scrub Sinks

SS Series Surgical Scrub Sinks

Combining the best materials, components, and craftsmanship available, each sink is hand built to the customer's specific need.

Standard Features:

- Deep sloping basin.
- Adjustable thermostatic controlled pressure regulating mixing valve, one per station.
- Check stop strainers for ease of maintenance.
- High mount swivel gooseneck spout with adjustable swivel spray aerator.
- Flow control valve.
- Backflow preventers.
- Infrared sinks equipped with 24 volt plug in transformer.
- Splash screen dividers with double and triple station sinks.
- Main support and Z-brackets (in-wall support optional).
- 304 stainless steel polished and ground to a #4 finish.



Knee Operated Surgical Scrub Sinks

The operation of these sinks is simple. Press the knee panel one time to activate water flow. Press again to stop flow. The soap panel works in a pumping fashion with soap being dispensed through a stainless steel spout located near the mixing valve. The base has been designed to accommodate gallon sized soap container storage behind the knee panels.



Infrared Operated Surgical Scrub Sinks

The infrared sensor is located on the front of each station and is activated by the user's waist. There is a 2 second on/off delay for water flow. The sensor range is 12" to 14". Both features combined eliminate activation by simply passing by.

Knee Operated Surgical Scrub Sink

Model# SS32-PB Pedestal Base Shown

Infrared Operated Surgical Scrub Sink

Model# SS32-IR-FSP



Deluxe Surgical Scrub Sink

Model# SS64-IR-MB

Deluxe Surgical Scrub Sinks

Built as a hybrid, these sinks have dual controls. Equipped with infrared operated water control as well as a manual bypass for knee operated water control. The main operation is via the infrared sensor located between the knee panels. The knee operation for water control is on standby as a backup. In addition the sink is also equipped with knee operated soap.

Tables, Sinks, and Warming Cabinets

Processing Sinks

Standard Features

- Fully-welded construction.
- Sound-proofed reinforced tops.
- Multiple layouts for individual requirements.
- All sinks are manufactured to individual specifications.
- Full CAD rendering of specified sink.



Sample Processing Sinks



Options

- Hydraulic adjustable heights for ergonomics.
- Full base unit underneath with door, drawer or open shelf combination.
- Pipe stand base unit available.
- Multiple bowls in any combination.
- Multiple faucet options including pre-rinse.
- Air hose with 5 blowgun tips.
- DI water.
- Foot operated drains.

Processing Sinks and Upper Cabinets Made to Order



Tables, Sinks, and Warming Cabinets

Blanket & Fluid Warming Cabinets

LogiQuip offers a variety of Blanket and Fluid Warming Cabinets. As with all of our products, warming cabinets undergo stringent testing procedures during the manufacturing process to assure a long life, the highest quality, and accurate temperature settings. We offer standard sizes along with custom sizes to meet specific requirements.

Standard Features

- Keyed temperature lock out.
- Push button digital operation.
- Large dual digital displays.
- Displays actual temperature and set point at all times.
- Separate digitally controlled heating chambers.
- Rapid warm time (2 - 6 hours depending on load).
- Uniform heating.
- Fully insulated.
- Audible and visual over-temp alarm.
- UL/CUL approved.
- FDA registered.
- OSHPD pre-approved.
- CSA Approved in Canada

Options

- Full view glass door.
- Mobile dolly or cart.
- 18" depth.
- Full recess.
- Pass through chamber.
- Left hand door hinge.
- Keyed door lock.
- Temperature recorder.
- Seismic brackets.
- 220/240 V, 50/60HZ.
- Custom sizes.



Model# TWC24 Overall Size

- 30"W x 26.5"D x 78.25"H.
- 14.00 cubic feet of storage.
- 1 adjustable shelf.

TWC24-G-MB
(shown on optional dolly)

Upper Chamber I.D.

- 26"W x 23"D x 11.25"H.
- 3.9 cubic feet of storage.
- Temperature range: 90-110 F.
- Accuracy +/- 1 degree F. of set point

Intermediate Chamber

- 26"W x 23"D x 10.75"H.
- 3.7 cubic feet of storage.
- Temperature range: 90-135 F.
- Accuracy +/- 1 degree F. of set point.
- Lower Chamber I.D.
- 26"W x 23"D x 18.50"H.

- 6.4 cubic feet of storage.
- Temperature range: 90-160 F.
- Accuracy +/- 1 degree F. of set point.

Electrical

- 120V, 13.5 amp., 50/60HZ.
- .65 kw/hr, 2238 btu/hr.
- UL/CUL listed.
- FDA registered.



Model# DWC24 Overall Size

- 30"W x 26.5"D x 74.5"H.
- 17.2 cubic feet of storage.
- 3 adjustable shelves.

DWC24-MB
(shown on optional dolly)

Upper Chamber I.D.

- 26"W x 23"D x 15.25"H.
- 5.27 cubic feet of storage.
- Temperature range: 90-160 F.
- Accuracy +/- 1 degree F. of set point.

Lower Chamber I.D.

- 26"W x 23"D x 34.5"H.
- 11.93 cubic feet of storage.
- Temperature range: 90-160 F.
- Accuracy +/- 3 degree F. of set point.

Electrical

- 120V, 12 amp., 50/60HZ.
- .75 kw/hr, 2557 btu/hr.
- UL/CUL listed.
- FDA registered.

Tables, Sinks, and Warming Cabinets

Blanket & Fluid Warming Cabinets



DWC2424-MB
(shown on optional dolly)

Model# DWC2424
Overall Size

- 24"W x 26.5"D x 74.5"H.
- 15.95 cubic feet of storage.
- 3 adjustable shelves.

Upper Chamber I.D.

- 20"W x 23"D x 15.25"H.
- 4.05 cubic feet of storage.
- Temperature range: 90-160 F.
- Accuracy +/- 1 degree F. of set point.

Lower Chamber I.D.

- 26"W x 23"D x 34.5"H.
- 11.9 cubic feet of storage.
- Temperature range: 90-160 F.
- Accuracy +/- 3 degree F. of set point.

Electrical

- 120V, 10 amp., 50/60HZ.
- .46 kw/hr, 1570 btu/hr.
- UL/CUL listed.
- FDA registered.



Model# DWC2460
Overall Size

- 24"W x 20.5"D x 64.75"H.
- 7.45 cubic feet of storage.
- 3 adjustable shelves.

Upper Chamber I.D.

- 20"W x 17"D x 13.5"H.
- 2.65 cubic feet of storage.
- Temperature range: 90-160 F.
- Accuracy +/- 1 degree F. of set point.

Lower Chamber I.D.

- 20"W x 17"D x 24.5"H.
- 4.80 cubic feet of storage.
- Temperature range: 90-160 F.
- Accuracy +/- 2 degree F. of set point.

Electrical

- 120V, 10 amp., 50/60HZ.
- .49 kw/hr, 1678 btu/hr.
- UL/CUL listed.
- FDA registered.

DWC2460-G

Model# DWC36
Overall Size

- 30"W x 26.5"D x 74.5"H.
- 17.59 cubic feet of storage.
- 4 adjustable shelves.

Upper Chamber I.D.

- 26"W x 23"D x 25"H.
- 8.63 cubic feet of storage.
- Temperature range: 90-160 F.
- Accuracy +/- 2 degree F. of set point.

Lower Chamber I.D.

- 26"W x 23"D x 26"H.
- 8.96 cubic feet of storage.
- Temperature range: 90-160 F.
- Accuracy +/- 2 degree F. of set point.

Electrical

- 120V, 12 amp., 50/60HZ.
- .69 kw/hr, 2356 btu/hr.
- UL/CUL listed.
- FDA registered.

DWC36-G-MB
(shown on optional dolly)



Model# DWC60
Overall Size

- 30"W x 20.5"D x 64.75"H.
- 9.73 cubic feet of storage.
- 3 adjustable shelves.

Upper Chamber I.D.

- 26"W x 17"D x 13.5"H.
- 3.46 cubic feet of storage.
- Temperature range: 90-160 F.
- Accuracy +/- 1 degree F. of set point.

Lower Chamber I.D.

- 26"W x 17"D x 24.5"H.
- 6.27 cubic feet of storage.
- Temperature range: 90-160 F.
- Accuracy +/- 2 degree F. of set point.

Electrical

- 120V, 10 amp., 50/60HZ.
- .49 kw/hr, 1678 btu/hr.
- UL/CUL listed.
- FDA registered.

DWC60-G



Tables, Sinks, and Warming Cabinets

Blanket & Fluid Warming Cabinets



SWC1518-G

Model# SWC1518 Overall Size

- 18"W x 15.5"D x 22"H.

Chamber I.D.

- 14"W x 13"D x 13"H.
- 1.4 cubic feet of storage.
- Temperature range: 90-135 F.
- Accuracy +/- 1 degree F. of set point.

Electrical

- 120V, 3.5 amp., 50/60HZ.
- .39 kw/hr, 1332 btu/hr.
- UL/CUL listed.
- FDA registered.

Model# SWC1824 Overall Size

- 24"W x 20.5"D x 24.5"H.

Chamber I.D.

- 20"W x 17"D x 15.25"H.
- 3.0 cubic feet of storage.
- Temperature range: 90-160 F.
- Accuracy +/- 1 degree F. of set point.

Electrical

- 120V, 3.5 amp., 50/60HZ.
- .45 kw/hr, 1537 btu/hr.
- UL/CUL listed.
- FDA registered.



SWC1824-G

Model# SWC24 Overall Size

- 30"W x 26.5"D x 24.5"H.
- 1 adjustable shelf.

Chamber I.D.

- 26"W x 23"D x 15.25"H.
- 5.27 cubic feet of storage.
- Temperature range: 90-160 F.
- Accuracy +/- 1 degree F. of set point.

Electrical

- 120V, 7.5 amp., 50/60HZ.
- .30 kw/hr, 1024 btu/hr.
- UL/CUL listed.
- FDA registered.



SWC24-G

Model# SWC36 Overall Size

- 30"W x 26.5"D x 36"H.
- 2 adjustable shelves.

Chamber I.D.

- 26"W x 23"D x 25"H.
- 8.63 cubic feet of storage.
- Temperature range: 90-160 F.
- Accuracy +/- 2 F. degree of set point.

Electrical

- 120V, 7.5 amp., 50/60HZ.
- .35 kw/hr, 1195 btu/hr.
- UL/CUL listed.
- FDA registered.



SWC36-G

Model# SWC2460 Overall Size

- 24"W x 20.5"D x 64.75"H.
- 3 adjustable shelves.

Chamber I.D.

- 20"W x 17"D x 49.5"H.
- 9.73 cubic feet of storage.
- Temperature range: 90-160 F.
- Accuracy +/- 1 degree F. of set point.

Electrical

- 120V, 7.5 amp., 50/60HZ.
- .49 kw/hr, 1678 btu/hr.
- UL/CUL listed.
- FDA registered.



SWC2460-G-MB
(shown on optional dolly)

Tables, Sinks, and Warming Cabinets

Blanket & Fluid Warming Cabinets



SWC60-G-MB
(shown above on optional dolly)

Model# SWC60 Overall Size

- 30"W x 20.5"D x 64.5"H.
- 3 adjustable shelves.

Chamber I.D.

- 26"W x 17"D x 51"H.
- 13.05 cubic feet of storage.
- Temperature range: 90-160 F.
- Accuracy +/- 3 degree F. of set point.

Electrical

- 120V, 7.5 amp., 50/60HZ.
- .46 kw/hr, 1570 btu/hr.
- UL/CUL listed.
- FDA registered.

Model# SWC72 Overall Size

- 30"W x 26.5"D x 74.5"H.
- 2 adjustable shelves.
- 2 fixed shelves.

Chamber I.D.

- 26"W x 23"D x 61"H.
- 21 cubic feet of storage.
- Temperature range: 90-160 F.
- Accuracy +/- 3 degree F. of set point.

Electrical

- 120V, 7.5 amp., 50/60HZ.
- .53 kw/hr, 1810 btu/hr.
- UL/CUL listed.
- FDA registered.



SWC72-G-MB
(shown above on optional dolly)

Optional Mobile Stands



T0007

T0023

Optional Mobile Base

T0087: 33"W x 28"D x 6.75"H.



* If ordered with the warmer, add -MB to the model #.

** Note this will add 2" to the standard warmer height.

Lock Options



Standard keyed
temperature lockout



Optional door lock

Temperature Monitoring and Recording Device, for Easy Integrated Data Logging

- User friendly to any computer with Excel 97 - 2010
- Each chamber is controlled independently
- USB access port on faceplate to allow for easy downloading
- Record temperature data for safe keeping
- Multiple warming cabinets can be controlled and monitored with one USB stick
- Remote alarm contact
- Keyed temperature lock out



Tables, Sinks, and Warming Cabinets

Blanket & Fluid Warming Cabinets

How to Order Blanket & Fluid Warming Cabinets

18" DEPTH	GLASS DOORS	MOBILE BASE	MOBILE STAND	RECESSED UNIT	LH DOOR	KEYED DOOR LOCK	TEMPERATURE RECORDER	220/240 V, 50/60 HZ	PASS THROUGH CHAMBER
Opt Ext	Opt Ext	Opt Ext	Opt Ext	Opt Ext	Opt Ext	Opt Ext	Opt Ext	Opt Ext	Opt Ext
18	G	MB	MS	R	LH	DL	TR	220	P

For Example: Model # DWC24 with the 18" depth option, Glass Doors, Recessed unit, Hinges on the left, Keyed door lock, 220/240 V 50/60 HZ

Order Model # DWC24-18-G-R-LH-DL-220

Standard Models

Single Chamber Units

MODEL #	W x D x H (inches)	CAPACITY cubic feet
SWC1518	18 x 15.5 x 22	1.4
SWC1518-G	18 x 15.5 x 22	1.4
SWC1824	24 x 20.5 x 24.5	3.0
SWC1824-G	24 x 20.5 x 24.5	3.0
SWC24	30 x 26.5 x 24.5	5.27
SWC24-G	30 x 26.5 x 24.5	5.27
SWC2460	24 x 20.5 x 64.75	9.73
SWC2460-G	24 x 20.5 x 64.75	9.73
SWC36	30 x 26.5 x 36	8.63
SWC36-G	30 x 26.5 x 36	8.63
SWC60	30 x 20.5 x 64.75	13.05
SWC60-G	30 x 20.5 x 64.75	13.05
SWC72	30 x 26.5 x 74.5	21
SWC72-G	30 x 26.5 x 74.5	21

Double Chamber Units

MODEL #	W x D x H (inches)	CAPACITY cubic feet
DWC24	30 x 26.5 x 74.5	17.2
DWC24-G	30 x 26.5 x 74.5	17.2
DWC2424	24 x 26.5 x 74.5	15.95
DWC2424-G	24 x 26.5 x 74.5	15.95
DWC2460	24 x 20.5 x 64.75	7.45
DWC2460-G	24 x 20.5 x 64.75	7.45
DWC36	30 x 26.5 x 74.5	17.59
DWC36-G	30 x 26.5 x 74.5	17.59
DWC60	30 x 20.5 x 64.75	9.73
DWC60-G	30 x 20.5 x 64.75	9.73

Triple Chamber Units

MODEL #	W x D x H (inches)	UPPER CAPACITY cubic feet	MIDDLE CAPACITY cubic feet	BOTTOM CAPACITY cubic feet
TWC24	30 x 26.5 x 78.25	3.9	3.7	6.4
TWC24-G	30 x 26.5 x 78.25	3.9	3.7	6.4

Note: Due to the vast number of sizes we have available, not all of our models are shown. Please let us know your specific dimensions to get the model that will best suit your needs.

If we don't currently have the size you need, we will build a custom one for you!

Warranty: All warming cabinets manufactured by LogiQuip will carry a lifetime guarantee against product craftsmanship, 1-year labor and 2-year parts guarantee. The factory will service all units without cost to the buyer for 1-year from shipment. After the 1-year period, replacement of a defective part (labor) will be at buyer's expense. We will exchange all defective parts at no cost to the buyer for a period of 2-years from shipment. All defective parts must be returned within 30 days to ensure proper credit. An RMA from LogiQuip. must be obtained prior to items return.

Tables, Sinks, and Warming Cabinets

Height Adjustable Reprocessing Work Station

Overhead light fixture

Adjustable shelf

Double bay frame

Stainless steel worksurface over sealed core

Hanging drawers

5" caster set
2 swivel/2 brake

Task light with magnification (3/4x)

Flat panel monitor holder

Horizontal power rail

Hand crank

Bottom shelf

CPU holder



Tables, Sinks, and Warming Cabinets



ALE/R6030-LOGI
Motorized Adjustment

Manual or Motorized Height-Adjustable Reprocessing Worktables... An Ergonomic Solution

- Height-adjustable from 31.5" low, to 47.5".
- Steel base with epoxy powder coat finish. Standard color: Light gray.
- Motorized (low voltage) or manual adjustment — Motorized has a 3-position memory for adjustment by each technician.
- 16-gauge #304 stainless steel sheet overlaid on particle board, 2"D.
- 500 lb. capacity.
- One step height adjustment. All work station accessory components raise and lower with the table automatically.
- 5" Casters come standard.

Height-Adjustable Work Tables with Stainless Steel Top

MODEL #	D x W x H	DESCRIPTION
	(inches)	
ALC/R4830-LOGI	30 x 48 x 31.5 - 47.5	Hand crank
ALC/R6030-LOGI	30 x 60 x 31.5 - 47.5	Hand crank
ALC/R4836-LOGI	36 x 48 x 31.5 - 47.5	Hand crank
ALC/R6036-LOGI	36 x 60 x 31.5 - 47.5	Hand crank
ALE/R4830-LOGI	30 x 48 x 31.5 - 47.5	Motorized
ALE/R6030-LOGI	30 x 60 x 31.5 - 47.5	Motorized
ALE/R4836-LOGI	36 x 48 x 31.5 - 47.5	Motorized
ALE/R6036-LOGI	36 x 60 x 31.5 - 47.5	Motorized

Height Adjustable Reprocessing Worktables

Table Mounted Upright Frames

MODEL #	DIMENSIONS	DESCRIPTION
	(inches)	
NXF/RM4848-2	48 x 48	Frame, double bay, rear justified
NXF/RM6048-2	60 x 48	Frame, double bay, rear justified
NXF/RM7248-2	72 x 48	Frame, double bay, rear justified

Shelves

MODEL #	DIMENSIONS	DESCRIPTION
	(inches)	
NXS/A2415L	15 x 24	Shelf, adjustable, 100 lb cap
NXS/A3015L	15 x 30	Shelf, adjustable, 100 lb cap
NXS/A3615L	15 x 36	Shelf, adjustable, 100 lb cap
NXS/V2418L	18 x 24	Shelf, variable angle, 150 lb cap
NXS/V3018L	18 x 30	Shelf, variable angle, 150 lb cap
NXS/V3618L	18 x 36	Shelf, variable angle, 150 lb cap

Drawer Systems

MODEL #	D x W x H	DRAWERS
	(inches)	
WBHC075-1	19.25 x 16.75 x 5	(1) 2.75"H
WBHC150-1	19.25 x 16.75 x 8	(1) 5.75"H
WBHC150-2	19.25 x 16.75 x 8	(2) 2.75"H
WBHC300-1	19.25 x 16.75 x 14	(1) File
WBHC300-2	19.25 x 16.75 x 14	(2) 5.75"H
WBHC450-2	19.25 x 16.75 x 20	(1) 5.75"H / (1) File
WBHC450-3	19.25 x 16.75 x 20	(3) 5.75"H

Includes key lock, bracket kit, all-steel construction, standard light gray, full-extension slides.

Additional Accessories

MODEL #	DESCRIPTION
CPU/H	CPU holder
NXA/MS-KB	Monitor shelf & keyboard holder
NX/MP	Mouse tray
NXRMT/3030-NT	Marker board, 2 sided
LTWL 0304	Task light with magnification (3/4X)
QQNX/AMB-BKT/EXT	Frame mounting bracket for LTWL 0304
QQNX/AMB-SLV	Frame mounting sleeve for LTWL 0304
NXHPS/30-15	Power rail, 30", horizontal, 15 AMP

OR & Reprocessing Storage Systems

Par Wall® and Par Stor-III®

Par Wall®

LEAN & CLEAN STORAGE

**Wall-Mount or Freestanding...
Dozens of Customized Configurations**

*Keep Supplies
Organized and Accessible*



OR Sterile Storage- Wall Mount Units

Par Stor-III®

HIGH-DENSITY

**All Vertical Space is Used.
Organized and Compartmentalized**

Wire Minimizes Dust



Row of Supplies in Reprocessing



Par Stor®, Par Wall®, and Container Carts in OR Core



Organized Single Instrument and General Storage

OR & Reprocessing Storage Systems

Par Stor-III®



PAR STOR-III

Par Stor® Complete Standard Module: Complete with 2 frames, 4 cross bars, top wire basket, and from top: 1-3"H basket on sloping slides with safety stop with 2 long and 2 short dividers; 5-5"H baskets each with 1 long and 2 short dividers; 1-5"H drawer with 2 long and 2 short dividers; 1-8"H basket with 1 long and 1 short divider; All on full-extension slides. Bottom dust cover and 4-4"H casters, 2 with brakes. 25 model #LH3 label tags included.


MODEL #	DESCRIPTION	D X W X H (OD)
	Par Stor® Complete Unit	(inches)
PAR STOR-III	Par Stor® Complete Unit, Starter, 7 Baskets, 1 Drawer	24 x 20 x 75
PARST ADD-III	Par Stor® Complete Unit, Add-On, 7 Baskets, 1 Drawer	24 x 19 x 75

* Add-On Units have all above components except with 1 frame only.

** Units can be customized using components listed on page 6.

Par Stor® Peel Pouch Storage Module:

Complete with 2 frames, 4 cross bars, 7 - 5"H wire baskets on full extension slides, 36 short dividers, 10 long dividers, 168 label tags, bottom dust cover and 4"H casters, two with brakes. Model #PSIR-III




Color Code Label Options




LH3 - Purple




LH3 - Blue



LH3 - Green




LH3 - Yellow




LH3 - Orange



LH3 - Red



LH3 - Pink



LH3 - Clear



Par Stor Container Storage Module:

Complete with 2 frames, 4 cross bars, 7 wire shelves on slides. Bottom dust cover and 4 - 4" swivel casters, two with brakes. Model #PS-CSS1-III



OR & Reprocessing Storage Systems

Par Wall®



PARWALL36-W (group of 4)

**A Simple Solution for Storing Large,
Irregular Shaped Items with Ease**

**Maximum Storage
Minimum Footprint**



P3619B Panel Attached to LP3674FSG
Panel Using PWBRACKET3619
Baskets/hook optional

Par Wall® Wall-Mount Complete Units- 6 Baskets:

Complete with 1 panel, 1 bottom solid dust shelf, 2-18"D large mesh baskets, 1-18"D fine mesh basket, 3-12"D fine mesh baskets. Each basket complete with mounting brackets, 5 front/back dividers and full-length label holders.

MODEL #	DESCRIPTION	D x W x H
	Par Wall® Complete Wall-Mount Unit	(inches)
PARWALL18-W	Complete Wall-Mount Unit, 6 Baskets	18 x 19 x 61
PARWALL36-W	Complete Wall-Mount Unit, 6 Baskets	18 x 37 x 61

Panel Mount Bracket Kits for Frame:

Epoxy coated steel cross frames for mounting louvered panels to freestanding frames. Brackets and hardware (nuts, bolts) included. Panels not included.

****Designed for LogiQuip panels only.****

MODEL #	DESCRIPTION
PWBRACKET1819	Par Wall® Panel/Frame Bracket 18"W (set of 2) for P1819B panel
PWBRACKET1861	Par Wall® Panel/Frame Bracket 18"W (set of 4) for P1861B panel
PWBRACKET3619	Par Wall® Panel/Frame Bracket 36"W (set of 2) for P3619B panel
PWBRACKET3661	Par Wall® Panel/Frame Bracket 36"W (set of 4) for P3661B panel



PARWALL36-INST Bulk Instrument
Storage for large, odd shaped instruments.
Hooks ordered separately.

OR & Reprocessing Storage Systems

Par Wall® Complete Units

6 Baskets

Par Wall® Single-Sided Stationary and Mobile Complete Units:

Complete with 1 frame, 1 bottom dust shelf, 2-18"D large mesh baskets, 1-18"D fine mesh basket, 3-12"D fine mesh baskets. Each basket complete with mounting brackets, 5 front/back dividers and full-length label holders. Add-on units have frames with 1 upright only.

MODEL #	DESCRIPTION	D x W x H
	Par Wall® Complete Unit	(inches)
PARWALL18	Par Wall® Complete Unit, 6 Baskets	18 x 18 x 74
PARWALL18M	Par Wall® Complete Unit, 6 Baskets, Mobile	26 x 18 x 78
PARWALL18-A	Par Wall® Complete Unit, 6 Baskets, Add-On	18 x 18 x 74
PARWALL18-AM	Par Wall® Complete Unit, 6 Baskets, Add-On, Mobile	26 x 18 x 78
PARWALL36	Par Wall® Complete Unit, 6 Baskets	18 x 36 x 74
PARWALL36M	Par Wall® Complete Unit, 6 Baskets, Mobile	26 x 36 x 78
PARWALL36-A	Par Wall® Complete Unit, 6 Baskets, Add-On	18 x 36 x 74
PARWALL36-AM	Par Wall® Complete Unit, 6 Baskets, Add-On, Mobile	26 x 36 x 78



PARWALL36M

12 Baskets

Par Wall® Double-Sided Stationary and Mobile Complete Units:

Complete with 1 double-sided frame, 2 bottom dust shelves, 4-18"D large mesh baskets, 2-18"D fine mesh basket, 6-12"D fine mesh baskets. Each basket complete with mounting brackets, 5 front/back dividers and full-length label holders.

MODEL #	DESCRIPTION	D x W x H
	Par Wall® Complete Unit	(inches)
PARWALL36-DB	Par Wall® Complete Unit, Dbl Sided, 12 Baskets	36 x 18 x 74
PARWALL36-DBM	Par Wall® Complete Unit, Dbl Sided, 12 Baskets, Mobile	40 x 18 x 78



PARWALL36-DBM

7 Baskets

Par Wall® Single-Sided Stationary and Mobile Complete Units:

Complete with 1 frame, 1 bottom dust shelf, 2-18"D large mesh baskets, 1-18"D fine mesh basket, 4-12"D fine mesh baskets. Each basket complete with mounting brackets (top basket to have slanted brackets) 5 front/back dividers and full-length label holders. Add-on units have frames with 1 upright only.

MODEL #	DESCRIPTION	D x W x H
	Par Wall® Complete Unit	(inches)
PARWALL187	Par Wall® Complete Unit, 7 Baskets	18 x 18 x 74
PARWALL187M	Par Wall® Complete Unit, 7 Baskets, Mobile	26 x 18 x 78
PARWALL187-A	Par Wall® Complete Unit, 7 Baskets, Add-On	18 x 18 x 74
PARWALL187-AM	Par Wall® Complete Unit, 7 Baskets, Add-On, Mobile	26 x 18 x 78
PARWALL367	Par Wall® Complete Unit, 7 Baskets	18 x 36 x 74
PARWALL367M	Par Wall® Complete Unit, 7 Baskets, Mobile	26 x 36 x 78
PARWALL367-A	Par Wall® Complete Unit, 7 Baskets, Add-On	18 x 36 x 74
PARWALL367-AM	Par Wall® Complete Unit, 7 Baskets, Add-On, Mobile	26 x 36 x 78



PARWALL367

* "M" mobile units complete with mobile kit and hardware. 4 x 4" swivel casters, two with brakes.

OR & Reprocessing Storage Systems

Conventional Storage Carts

LogiQuip has a huge range of conventional wire products. Here is a glimpse of popular mobile wire shelving units used in OR and reprocessing. For our full range of products please refer to our full wire catalog.

Standard Unit Components

- 4 each 74"H posts.
- 5 wire shelves with split sleeves.
- 5" swivel non-marking polyurethane casters.
- Donut bumpers.

"Complete" Unit Accessories (-CL)

- 5 full-length 3-position label holders.
- Bottom shelf cover, solid, flat, poly.

Weight Capacity

- 5" casters. 250 lbs. capacity each.
- Total cart 750 lbs. capacity.

Other Standard Sizes Avail.

- Depths: 14" and 21"
- Widths: 24", 30", 42" and 54"
- Posts: 54", 63" and 86"

Open Units

5 Shelf


Complete with

4 Swivel Casters (2 with brakes),
1 solid bottom shelf cover and
5 full-length 3-position label holders.

D x W x H (inches)	WT. (lbs.)	CHROME PLUS Part #	DURO SEAL Part #
18 x 36 x 80	75	MB835CH-CL	MB835PW-CL
18 x 48 x 80	89	MB855CH-CL	MB855PW-CL
18 x 60 x 80	110	MB865CH-CL	MB865PW-CL
18 x 72 x 80	117	MB875CH-CL	MB875PW-CL
24 x 36 x 80	94	MB435CH-CL	MB435PW-CL
24 x 48 x 80	113	MB455CH-CL	MB455PW-CL
24 x 60 x 80	138	MB465CH-CL	MB465PW-CL
24 x 72 x 80	149	MB475CH-CL	MB475PW-CL



MB455CH-CL



MEB835CH-CL

Enclosed units

5 Shelf

Complete with

4 Swivel Casters (2 with brakes),
1 solid bottom shelf cover and
5 full-length, 3-position label holders.

D x W x H (inches)	WT. (lbs.)	CHROME PLUS Part #	DURO SEAL Part #
18 x 36 x 80	115	MEB835CH-CL	MEB835PW-CL
18 x 48 x 80	137	MEB855CH-CL	MEB855PW-CL
18 x 60 x 80	165	MEB865CH-CL	MEB865PW-CL
18 x 72 x 80	178	MEB875CH-CL	MEB875PW-CL
24 x 36 x 80	141	MEB435CH-CL	MEB435PW-CL
24 x 48 x 80	168	MEB455CH-CL	MEB455PW-CL
24 x 60 x 80	200	MEB465CH-CL	MEB465PW-CL
24 x 72 x 80	217	MEB475CH-CL	MEB475PW-CL

OR & Reprocessing Storage Systems

TSUN is our premium wire shelving line: Smooth waterfall edge on front of shelf eliminates ripping of wrapped goods. Shelves adjust without any tools. Built-in endscreens. TSUN is also compatible with most standard wire accessories (refer to our "Storage and Transport Product Guide")

A MUST FOR EVERY STERILE STORAGE AREA, ESPECIALLY INSTRUMENT SETS!

TSUN

Complete TSUN Cart Consists of:

Backscreen of
3" x 3" wire mesh

**READILY
ADJUSTABLE
NO TOOLS
REQUIRED!**

Chrome plated wire
shelves each with flat-face
label holders.

Metal Solid Bottom
Shelf Cover

TSUN2464-5H-CL

4 x 5" poly swivel casters,
two with brakes

3 joining bars (1 at
top back, 2 at bottom)



Two end frames function
as enclosure panels

TSUN Complete Carts

MODEL #	DESCRIPTION	D x W x H (inches)
TSUN2436-5SH-CL	Complete TSUN 5-Shelf Unit 24" x 36"	28 x 40 x 73
TSUN2448-5SH-CL	Complete TSUN 5-Shelf Unit 24" x 48"	28 x 52 x 73
TSUN2464-5SH-CL	Complete TSUN 5-Shelf Unit 24" x 64"	28 x 68 x 73
TSUN2472-5SH-CL	Complete TSUN 5-Shelf Unit 24" x 72"	28 x 76 x 73

MODEL #	D x W (inches)
TSUN2436-PL	TSUN Plastic Shelf Cover 24 x 36
TSUN2448-PL	TSUN Plastic Shelf Cover 24 x 48
TSUN2464-PL	TSUN Plastic Shelf Cover 24 x 64
TSUN2472-PL	TSUN Plastic Shelf Cover 24 x 72

MODEL #	DESCRIPTION (inches)
CFDIV246	Shelf Divider 6"H
CFDIV248	Shelf Divider 8"H
TSUN2436-SH	TSUN Wire Shelf 24" x 36"
TSUN2448-SH	TSUN Wire Shelf 24" x 48"
TSUN2464-SH	TSUN Wire Shelf 24" x 64"
TSUN2472-SH	TSUN Wire Shelf 24" x 72"

OR & Reprocessing Storage Systems

Suture/Lens/Implant Carts



SL536CH

Single Sided Suture/Lens/Implant Carts Include

- 4 posts 63"H.
- 4 (or 5- see table) slanted baskets each with label holders. **NOT PICTURED**
- Bottom wire shelf with solid shelf cover. (for infection control)
- 4-5" swivel casters with donut bumpers, two with brakes.

MODEL #	DESCRIPTION	D x W x H (OD)
	Single Sided	(inches)
SL424CH	4-Level 24" Single-Sided Suture Cart	18 x 24 x 69
SL436CH	4-Level 36" Single-Sided Suture Cart	18 x 36 x 69
SL448CH	4-Level 48" Single-Sided Suture Cart	18 x 48 x 69
SL524CH	5-Level 24" Single-Sided Suture Cart	18 x 24 x 69
SL536CH	5-Level 36" Single-Sided Suture Cart	18 x 36 x 69
SL548CH	5-Level 48" Single-Sided Suture Cart	18 x 48 x 69
SL24CH	Individual 24" Single Sided Suture Basket	18 x 24 x 10
SL36CH	Individual 36" Single Sided Suture Basket	18 x 36 x 10
SL48CH	Individual 48" Single Sided Suture Basket	18 x 48 x 10



STR24CH

Double Sided Suture/Lens/Implant Carts Include

- 4 posts 63"H.
- 4 (or 5- see table) slanted baskets each with label holders on both sides. **NOT PICTURED**
- Bottom wire shelf with solid shelf cover. (for infection control). Top wire shelf.
- End enclosure panels of 3" x 3" wire mesh
- 4-5" swivel casters with donut bumpers, two with brakes.

MODEL #	DESCRIPTION	D x W x H (OD)
	Double Sided	(inches)
STR424CH	4-Level 24" Double-Sided Suture Cart	24 x 24 x 69
STR436CH	4-Level 36" Double-Sided Suture Cart	24 x 36 x 69
STR448CH	4-Level 48" Double-Sided Suture Cart	24 x 48 x 69
STR524CH	5-Level 24" Double-Sided Suture Cart	24 x 24 x 69
STR536CH	5-Level 36" Double-Sided Suture Cart	24 x 36 x 69
STR548CH	5-Level 48" Double-Sided Suture Cart	24 x 48 x 69
STR24CH	Individual 24" Double Sided Suture Basket	24 x 24 x 69
STR36CH	Individual 36" Double Sided Suture Basket	24 x 36 x 10
STR48CH	Individual 48" Double Sided Suture Basket	24 x 48 x 10



STR524CH


OR & Reprocessing Storage Systems

Basket/IV Solution Carts- Highest Density in the Market

- 4 posts 74"H.
- Open front wire baskets & shelves. "SB style".
- Each basket complete with label holders and set # of dividers.
- Bottom wire shelf with solid shelf cover. Top wire shelf.
- 4-5" swivel casters with donut bumpers, two with brakes.



DV2010CH



CFDIV246



AF411CH



AF020CH



AF Basket Cart System

18" Deep		24" Deep		
MODEL #	D x W x H (OD) (inches)	MODEL #	D x W x H (OD) (inches)	DESCRIPTION
				Complete IV Cart
AF008CH	18 x 24 x 80	AF2008CH	24 x 24 x 80	Complete IV Cart 24"- 2 shelves, 6 baskets, 12 dividers
AF009CH	18 x 36 x 80	AF2009CH	24 x 36 x 80	Complete IV Cart 36"- 2 shelves, 6 baskets, 18 dividers
AF010CH	18 x 48 x 80	AF2010CH	24 x 48 x 80	Complete IV Cart 48"- 2 shelves, 6 baskets, 24 dividers
DV2010CH	18 x 6	N/A	N/A	Wire Basket Divider Standard, 18"
N/A	N/A	CFDIV246	24 x 6	Wire Basket Divider Standard, 24"
N/A	N/A	CFDIV248	24 x 8	Wire Basket Divider Standard, 24"

18" Deep		24" Deep		
MODEL #	D x W x H (OD) (inches)	MODEL #	D x W x H (OD) (inches)	DESCRIPTION
				Combo Cart
AF408CH	18 x 24 x 80	AF2408CH	24 x 24 x 80	Combo Cart 24"- 3 shelves, 4 baskets, 8 basket dividers
AF508CH	18 x 24 x 80	AF2508CH	24 x 24 x 80	Combo Cart 24"- 2 shelves, 5 baskets, 10 basket dividers
AF409CH	18 x 36 x 80	AF2409CH	24 x 36 x 80	Combo Cart 36"- 3 shelves, 4 baskets, 12 basket dividers
AF509CH	18 x 36 x 80	AF2509CH	24 x 36 x 80	Combo Cart 36"- 2 shelves, 5 baskets, 15 basket dividers
AF410CH	18 x 48 x 80	AF2410CH	24 x 48 x 80	Combo Cart 48"- 3 shelves, 4 baskets, 16 basket dividers
AF510CH	18 x 48 x 80	AF2510CH	24 x 48 x 80	Combo Cart 48"- 2 shelves, 5 baskets, 20 basket dividers
AF020CH	18 x 48 x 80	AF2020CH	24 x 48 x 80	Combo Cart 48"- 2 shelves, 5 baskets, 10 basket dividers
AF411CH	18 x 48 x 80	AF2411CH	24 x 48 x 80	Combo Cart 48"- 3 shelves, 4 baskets, 20 basket dividers
AF511CH	18 x 48 x 80	AF2511CH	24 x 48 x 80	Combo Cart 48"- 2 shelves, 5 baskets, 25 basket dividers


OR & Reprocessing Storage Systems

Floor Track System

Floor Track System Benefits


- CAD layout options at no cost or obligation.
- We will not recommend a system unless users gain at least 20% more linear footage.
- Units custom built to fit all shelving, existing or new. (Save \$\$\$ using existing)
- Choose from light-duty skates for individual unit tracks or heavy-duty carriages for multiple, back-to-back or mechanical assist units.
- No subfloors construction needed.
- Completely installed, including track install, transfer and assembly of existing units.
- System can easily be added-to and transferred with most parts re-used.

BEFORE



Traditional Static Storage
420 Lineal Feet Capacity

AFTER



Active Aisle System, 600 Linear Feet,
Increased Capacity 43% Over Static System

OR & Reprocessing Storage Systems

Floor Track System

Carriage System

- Carriages can be up to 12' in length and 4' in depth.
- Ball bearing design allows units up to 2000 lbs. (907kg) to move easily.
- "V" groove floor track is low-profile, cart-friendly and non-trip.
- Carriages: Satin anodized aluminum rails with roller bearing wheel assemblies.
- Track: Aluminum-extruded "V"-design (3/8" at highest peak) with hardened steel inserts. Reduces friction and provides positive carriage retention.
- Low profile end stops.
- Optional Mechanical Assist: Ergonomically designed easy move aisles, complete with built in aisle control.



Skate System

- Simple and cost-effective solution.
- Easy to configure: Determine length of track needed and provide depth of carts used.
- Aluminum skates with HDP end stops and roller bearing wheels attach to bottom of posts.
- "V" groove floor track is low-profile, cart-friendly and non-trip. (same style as carriages) 3/8"H.



Shown with User-Friendly Mechanical Assist



Low-Cost High-Density Storage on Skates

OR & Reprocessing Storage Systems

LogiCell® Room Carts

Fully Loaded for each OR Service


- Carts come complete with pictured accessories.
- Dividers included with most tubs, baskets and shelves.
- Units shown with tambour door.
- Hinged glass doors available on double & triple carts.
- Double carts: 28"D x 40"W x 80"H.
- Triple carts: 28"D x 57"W x 80"H.



2652VASC
Vascular Cart



2652CYSTO
Cysto Cart



2652CARD
Cardiac Cart



2652ARTH
Arthroscopy Cart




2652NEU
Neuro Cart



2652OPHENT
Ophthalmology/ENT Cart



2652LAP
Laparoscopy Cart



2652ORTH
Ortho Cart

OR & Reprocessing Storage Systems

LogiCell® Room Carts



2653VASC
Vascular Cart



2653ARTH
Arthroscopy Cart



2653CYSTO
Cysto Cart



2653CARD
Cardiac Cart



2653PACK
Drape/Pack Cart



2653STAP
Staple Cart



2653ORTH2
Ortho Cart



2653ORTH1
Ortho Cart



2653LAP
Laparoscopy Cart

OR & Reprocessing Storage Systems

QT Carts /Equipment Carts

Standard QT Carts

Ready for Quick Ship



34QT Series

20"D x 27"W x 34"H overall
2-3"H, 1-6"H, 1-9"H drawer
Drawer size: 15.75"D x 24"W



37QTN Series

21"D x 20"W x 37"H overall
2-3"H, 2-6"H, 1-9"H drawer
Drawer size: 15.75"D x 14"W



40QT Series

20"D x 27"W x 40"H overall
2-3"H, 2-6"H, 1-9"H drawer
Drawer size: 15.75"D x 24"W

Equipment Cart

Ready for Quick Ship



12" Adjustable Range

PK3366-WHT

24"D x 20"W x 34.5" - 46.5"H
2-3"H, 2-6"H drawer
Drawer size: 10.75"D x 14"W
Worktop size: 14"D x 18"W

Front Adjustable Pedal

Modular Divider Sets

One 3" H and one 6" H included with each cart



Model # QT-AD3

3" Adjustable divide set

Model # QTN-AD3

Narrow cart set

(3 x short + 6 x adjustable + 2 x 3" strip
+ 6 x T clips)



Model # QT-AD6

6" Adjustable divide set

Model # QTN-AD6

Narrow cart set

(3 x short + 4 x adjustable + 2 x 6" strip
+ 4 x T clips)

OR & Reprocessing Storage Systems

QT Carts /Equipment Carts

QT Cart Special Features

Included with each Cart



Soft Return, silent closing drawers



ABS top surface
Light, easy to clean

Sliding table top extension
Max load 17 lbs....



5" Medical grade quiet poly caster w/brake



Central lock
secures all drawers



Optional eLock System

Stock Accessories (Ready for immediate shipment)



VMS *
Multi storage
basket and rail



VFH *
File holder and rail



VGD *
Glove box holder and rail



VIV
IV pole



CPD
Computer drawer




P-AB Trellis
Trellis



VBS-1
Small fixed bins



VTB-1
Large tilt bins
(5 Bins)



VTB-2
Medium tilt bins
(12 Bins)



VGB *
Garbage can w/ Lid
VWB *
Open waste bin
(not pictured)

*** Accessories Requiring #VAB Accessory Bar
VAB Bar included with item**

Terms & Conditions

Acceptance

All orders subject to acceptance of Terms and Conditions of sale, regardless of the form in which submitted. All orders are quoted F.O.B. Galesburg, Michigan unless stated otherwise. Custom manufactured products may be subject to a deposit at the time of order.

Order Confirmation

To avoid duplication of shipment, the word "CONFIRMATION" or "CONFIRMING" should appear on all documents which are sent to confirm orders placed by telephone, fax or any other way.

Terms

Minimum order \$50.00. Terms are Net 30 days (subject to credit approval). An annual service charge of 24% will be assessed on past due accounts.

Shipments

Shipments will be made F.O.B. shipping point and the buyer assumes responsibility for all damages in transit. Such damages must be declared at time of receipt, directly on the shipping documents, and claims must be filed immediately with the transport company. If after receipt of shipment, concealed damage is discovered, immediately notify the delivering carrier requesting an inspection. The inspection report should be used to support your claim against the carrier. Any claim for damage or shortage must be filed by you with the delivering carrier. Unless a preferred routing is specified at the time of order, shipment will be made the best way, transportation prepaid and added to invoice.

Returned Goods Policy

1. All returns must be authorized.
2. Each request must include the following:
 - A. Customer name and address.
 - B. Invoice number
 - C. Invoice date
 - D. Customer P.O. number
 - E. Quantity, catalog number and description of item
 - F. Reason for return
3. Upon receipt of your request, a return authorization number will be issued and this must be clearly referenced with your return.
4. The return shipment must be coordinated with a customer service representative, who will advise the proper mode of transport for the return.

5. Returned goods may be subject to a minimum 20% handling charge based on state, condition of goods and reason for return.
6. Returns will NOT be authorized for the following:
 - A. Items purchased on a "special order" basis (unless original vendor will accept the return).
 - B. Items that are custom manufactured.
 - C. Partial units of sale.
 - D. Items without original packaging materials.
7. RETURNS MUST BE PREPAID AND WITHIN 10 DAYS OF ORIGINAL SHIPMENT.

Specifications

All dimensions and shipping weights are approximate. LogiQuip reserves the right to change specifications or substitute products of equal specifications without notice. The company also reserves the right to cease manufacturing or distribution of any product without notice. Manufacturer's specifications and prices are subject to change without notice.

Freight Classification

Freight Class is product specific and may vary from CL70 to CL250.

Installation

Installation of cabinets: Must be in compliance with all state and local code regulations and should be performed by experienced and qualified tradespeople. For best long term results and proper functioning of the cabinet doors and latches, the cabinets must be installed plumb and squarely true with adjoining cabinets. Compliance with all state and local regulations relative to anchoring cabinets should be confirmed by the appropriate hospital staff. Improper installation, or failure to comply with state and local regulations, may result in failure of the product and denial of potential warranty claims.

Warranty Statement

LogiQuip equipment is warranted to be free from defective material and workmanship under normal use and service for a period of one year from date of shipment to customer. During this period, LogiQuip will repair or replace (at its option) without charge, F.O.B. plant of manufacture, the equipment or any part thereof which factory inspection discloses to be defective. Shipments of defective material for repair must be sent prepaid.

Index

A

Accessories, Case Carts	5, 8
Adjustable Work Tables	36-37
Assembly Carts.....	11

B

Back Tables	13
Basket/IV Carts	45
Basin Stands.....	16
Blanket Warming Cabinets.....	31-35
Buckets	16

C

Carts, QT Specialty.....	50-51
Cabinets, Mobile, SS	22
Cabinets, Nurse Doc Station	26
Cabinets, Scope	27
Cabinets, SS Storage	22-27
Cabinets, Stationary.....	23, 24
Cabinets, Wall & Base	25
Cabinets, Warming.....	31-35
Case Carts, Closed	2-5
Case Carts, Open.....	6-8
Case Carts, O-TR	9
Casters, SS Case Cart.....	5
Containers, Wash Carts.....	10

D

Decontamination Sinks.....	30
----------------------------	----

F

Fluid Warming Cabinets	31-35
Floor Track Storage System	46, 47

H

Hampers, Linen.....	20
High Density Storage System, Floor Track	46-47

I

Instrument-Back Tables	13
IV Poles	19

K

Kick Buckets	16
--------------------	----

L

Label Holders, Par Stor	39
Label Holders, Par Wall.....	39
Label Holders, Wire	39
Label Tags, Colored	39
Lens Carts	44
Linen Hampers.....	20
LogiCell Carts.....	48-49

M

Mayo Stands.....	17
------------------	----

N

Nesting Back Tables	15
Nurse Documentation	26

O

Ortho/Implant Carts.....	12, 44, 48, 49
Over-The-Road/Offsite Carts	9

P

Processing Sinks	30
Par Cart System	12
Par Stor Storage System	39
Par Wall Storage System	40, 41
* Wall-Mount Complete Systems.....	41
Par Stor Modules	39
Peel Pouch Storage	39
Pick Carts	11
Platforms.....	18

R

Ring Stands.....	16
------------------	----

S

Scope Cabinets	27
Scrub Sinks	28, 29
Seating, Surgeon	18
Shelving, Mobile Wire	11, 42-45
Shelving, Wire	42-45
Sinks, Processing	30
Slanted Shelf Carts.....	44
Solution/IV Carts	45
Step Stools.....	18
Suture Carts	44

T

Tables, Back	13
Tables, Utility	14
Tables, Work	15, 36, 37
Tables, Adjustable	36, 37
Track System	46, 47
TSUN Wire System	43

U

Utility Tables.....	14
---------------------	----

W

Warmers, Blanket/Fluid.....	31-35
Wash Carts, Containers.....	10
Wall-Mount Supply Storage, Par Wall	40
Work Table (Stations).....	15, (36, 37)
Wire Shelves.....	11, 42-47
Wrap Racks	11



Healthcare Storage Solutions

5466B Timberlea Blvd.
Mississauga, Ontario, Canada L4W 2T7
1-800-551-0560 Fax: 888-625-0151

