



# LOGIQUIP

Healthcare Storage Solutions

**Anesthesia Carts**

**Emergency Carts**

**Isolation Carts**

**Bedside Carts**

**Pediatric Carts**

**Specialty Carts**

**Treatment Carts**



**QT**® **SERIES**

**THE QUIET TREATMENT CARTS**


*Quiet Design* with Your Patient in Mind

# QT®Cart Special Features

Included with each Cart



Soft Return, silent closing drawers



ABS top surface  
Light, easy to clean



Sliding table top extension  
Max load 17 lbs.

5" Medical grade quiet poly caster w/brake



## QUIET FEATURES

- Soft return, silent closing drawers
- Silent, hospital grade casters
- Quiet, innovative eLock System

## UNIQUE STANDARD FEATURES

- Modern clean design
- Larger capacity drawers in standard footprint
- Simple accessory rail system
- Lightweight aluminum frame
- Integrated bumper system
- (1) 3" & (1) 6" modular divider set
- ABS top with raised lip
- Pull out side shelf

## IDEAL FOR:

**Patient Care Areas:** ICU, CCU, Pediatric, Pre Op, PACU, Isolation and more

**Emergency Dept:** Laceration, ENT, Suture, General Treatment, Trauma and Airway

**Oncology:** Chemo and Infusion


**Phlebotomy:** Blood Draw

**Labs:** Endoscopy, Radiology, and Cath Lab

### 41QT Series

Shown with accessories in standard white

**Drawer Front Color Options** (Stocked Standard in white)



Central lock secures all drawers



Optional eLock System with proximity switch

**LOCK OPTIONS**

- Central key lock – standard
- Breakaway lock
- Keyless electronic lock
- HID Prox lock
- HID iCLASS lock

**ELOCK SPECS.**

- Whisper quiet operation
- 3,000 user capability
- Automatic relock
- Low battery indicator
- Key override
- Dual credential capability
- Available Wi-Fi allows IT to remotely add/delete users, automatically download audit trails, monitor battery level, and more



**VLS**  
Lock Management Software



**VPT**  
HID Proximity Tag



**VPC**  
HID Proximity Card

# Standard QT® Carts

\*Ready for Quick Ship



**QT35W21K3369W**

**\*35QT Series**

**Overall Dimensions:**  
20"D x 27"W x 34.5"H

**Drawers:**  
2-3", 1-6", 1-9"

**Drawer size:**  
15.75"D x 24"W



**QT38W24K3669W**

**38QT Series**

**Overall Dimensions:**  
20"D x 27"W x 37.5"H

**Drawers:**  
1-3", 2-6", 1-9"

**Drawer size:**  
15.75"D x 24"W



**QT41W27K33669W**

**\*41QT Series**

**Overall Dimensions:**  
20"D x 27"W x 40.5"H

**Drawers:**  
2-3", 2-6", 1-9" drawer

**Drawer size:**  
15.75"D x 24"W



**QT44W30K333669W**

**\*44QT Series**

**Overall Dimensions:**  
20"D x 27"W x 43.5"H

**Drawers:**  
3-3", 2-6", 1-9" drawer

**Drawer size:**  
15.75"D x 24"W

*Perfect for bedside use*



**QTD34W21K3666W**

**34QTD Series**

**Overall Dimensions:**  
21.5"D x 44"W x 34"H

**Drawers:**  
1-3", 3-6"

**Drawer Size:**  
Left: 15.75"D x 24"W  
Right: 15.75"D x 14"W



**QTN40W27K33669W**

**\*40QTN Series**

**Overall Dimensions:**  
21"D x 20"W x 39.5"H

**Drawers:**  
2-3", 2-6", 1-9"

**Drawer Size:**  
15.75"D x 14"W



## Specialty Treatment Carts and Accessory Packages



**QT41R27B33669R**

Shown with Accessories

### QT-EMERG - PKG

- VIV - IV Pole
- VDS - Defibrillator Shelf
- VCPR - Cardiac Board and Hooks
- VO2 - Oxygen Tank Holder
- VTR - Three-Sided Rail
- AD3(x2) - Adjustable Drawer Dividers, 3"
- AD6 - Adjustable Drawer Dividers, 6"
- TOS - 15 - Outlet Strip



**QT44W30BPED52**

### QT-ANES - PKG

- VIV - IV POLE
- P-AB-Back Bridge
- (2) VTB-1 - Tilt Bins, 5-Bin System
- VGD - Glove Dispenser
- VWB - Waste Basket
- AD3 - Adjustable Drawer Dividers, 3"
- AD6 - Adjustable Drawer Dividers, 6"
- VTR - Three-Sided Rail
- TD-1 - Tape Dispenser



**QT44Y30K3999Y**

Shown with Accessories

### QT-ISO - PKG

- VWB - Waste Basket
- VGB - Glove Dispenser
- VTR - Three-Sided Rail
- VUH - Utility Hooks
- VCR - Coat Rack



**QT41B27K33669B**

Shown with Accessories

**Stock Accessories** (Ready for immediate shipment)



**VAD-1**  
Locking Med Drawer



**CPD**  
Computer drawer

**Modular Divider Sets**

One 3" H and one 6" H included with each standard cart




**Model # QT-AD3**  
3" Adjustable divider set  
**Model # QTN-AD3**  
Narrow cart set  
(3 x short + 6 x adjustable  
+ 2 x 3" strip + 6 x T clips)



**Model # QT-AD6**  
6" Adjustable divider set  
**Model # QTN-AD6**  
Narrow cart set  
(3 x short + 4 x adjust-  
able + 2 x 6" strip  
+ 4 x T clips)



**Model VMT-1**  
Med Tray  
15"x22"x2"  
9 Compartments



**VTB-2 ‡**  
Medium tilt bins  
(12 Bins)

**VTB-1 ‡**  
Large tilt bins  
(5 Bins)

## Stock Accessories Continued (Ready for immediate shipment)

‡ Requires Trellis




**VO2**  
O<sup>2</sup> Holder



**VUH**  
Utility Hook Rack



**VIV**  
IV pole



**VBS-1 ‡**  
Small fixed bins  
each bin: 7"x4"x3"



**VCPR**  
Cardiac Board



**VDS**  
Defibrillator Shelf



**VGB \***  
Garbage can w/ Lid



**VSC**  
Sharps Container Set



**TD-1**  
Tape Holder



**VCR**  
Coat Rack



**P-AB Trellis**  
Trellis



**VTR**  
Three Sided Rail



**VGD3**  
Triple Glove  
Box Holder and rail



**VFH**  
File holder and rail  
6"x15"x22"




**VGD**  
Glove box holder and rail



**VMS**  
Multi storage basket & rail  
7"x13"x5"





Check Out Our Full Line of Mobile Healthcare Storage Solutions



 **LOGIQUIP**  
Healthcare Storage Solutions

P.O. Box 278  
Galesburg, Michigan 49053  
[www.logiquip.com](http://www.logiquip.com) [info@logiquip.net](mailto:info@logiquip.net)  
Phone: 800-665-3760  
Fax: 888-965-3984