



LOGICELL[®]

A World of



Modular Possibilities



LOGIQUIP
Healthcare Storage Solutions

About Us

Our History, Mission and People



Our Story

LogiQuip has been in the healthcare storage solutions business since 1992 with principals that have been experts in the field for a combined 90 years.

Historically, our main departments were OR/SPD and Materials Management. The new LogiCell[®] line brings us into the exciting Cath Lab and Interventional Radiology areas.

LogiQuip has partnered with FlairWood Industries, Inc. to produce the unique patented cell system design in carts and cabinets. From the marriage of experience and innovation we give to you the LogiCell[®] System of products.

Our Mission

Because each hospital's inventory, supply, and storage requirements are similar in theory, yet diverse in practice and application, our mission is to provide true solutions for inventory, space, storage, and material handling problems on an individual basis. The solution to a specific problem may consist of a single cart or a complete supply storage system! Whatever the answer, LogiQuip makes it workable and efficient.


Our mission: Solving storage and distribution issues with the best products and designs possible. We believe "space" to be one of the hospital's most expensive capital assets; as such, our philosophy is to use every cubic inch of that very expensive space.

We focus on maintaining infection control standards and functional ergonomics, which are key areas in healthcare storage and distribution.


Our products work in synergy: What we've designed for the OR can work in the Cath Lab and vice versa. We have the most varied line of storage and distribution products for all areas of the hospital, from the receiving dock to the OR.

Our People

We have in-depth experienced representation across North America. A smiling storage solutions adviser can meet with you at your convenience. We have highly trained customer service and CAD layout specialists in our Michigan home office. No job is too big or too small.



Typical SPD Layout Before: 1600 linear feet for storage



LogiQuip SPD Layout After: 2700 linear feet for storage

Table of Contents


A Directory to Your Modular World

The LogiCell [®] Advantage	2
LogiCell [®] Supply Carts	3 - 10
LogiCell [®] Tambour Door Cabinetry	11 - 15
LogiCell [®] Casework, Laminate & Stainless	16 - 21
LogiCell [®] Scope Cabinets	22 - 23
LogiCell [®] & Other Accessories	24 - 26
Wire Products	27 - 29
Samples, In Use	30 - 32
Terms and Conditions	33



The LogiCell® Advantage

What is the LogiCell® Advantage?



Innovative Side Panels

- ABS plastic with integrated safety stop
- Full access, 80% extension
- Slant standard baskets, tubs and shelves. No additional parts required
- Built into carts and cabinets, not bolted on
- Unique cell system is patent applied for




Unique Tambour Door System

- Designed with internal structural supports
- Cassette installation system allows for easy repair and replacement



Clean and Sealed Hinged Doors

- On stainless cabinetry and double tall carts
- Aluminum, not wood, with safety glass
- Flush-mount with no gaps



High Density Baskets and Tubs

- Long and short dividers in wire and plastic
- Baskets dividable in 1" increments
- Tubs dividable in 4" increments



Cabinetry Choices

Wide variety of options in laminate or stainless steel. Many configurations are achievable.



LogiCell® Supply Carts

LogiCell® Cart Family

Endless Supply Possibilities with LogiCell® Carts

- Constructed in 20-gauge textured epoxy coated steel with Frosty White finish
- Cell system of plastic side panels (patent pending) accommodates a range of accessories
- Standard tambour door with cassette design
- Tall double carts available with hinged glass doors of durable aluminum with hospital-grade safety glass
- Tubular handle both ends. Wrap-around bumper. 5" swivel casters, two with brakes



Tall Double Cart

- Overall: 28"D x 40"W x 80"H
- Tambour or hinged glass door

Tall Triple Cart

- Overall: 28"D x 57"W x 80"H
- Tambour door only

Tall 3/4 Cart

- Overall: 28"D x 40"W x 63.75"H
- Tambour or hinged glass door


Half Double Cart

- Overall: 28"D x 40"W x 42"H
- Tambour door only




LogiCell® Supply Carts

Half Cart and 3/4 Cart Shell Configurations



LogiCell® Supply Carts

Double Cart Shell Configurations



Model #2652
Tall Double Cart with 2 columns



Electronic Lock

Model #2651SH
Tall Double Cart with shelves



Model #2652TS
Tall Double Cart with Top Split



Model #2652BS
Tall Double Cart with Bottom Split



LogiCell® Supply Carts

Pre-Configured Double Carts for OR

Fully Loaded for each OR Service

- Carts come complete with pictured accessories
- Dividers included with most tubs, baskets and shelves
- Customize your own cart using shells from pages 4, 5 and 7 and accessories from pages 24 - 26
- Units shown with tambour door. Hinged glass doors available on double carts.
- Double carts: 28"D x 40"W x 80"H
- Triple carts: 28"D x 57"W x 80"H




Model #
2652VASC
Vascular Cart



Model #
2652CYSTO
Cysto Cart



Model #
2652CARD
Cardiac Cart



Model #
2652ARTH
Arthroscopy Cart



Model #
2652NEU
Neuro Cart



Model #
2652OPHENT
Ophthalmology/ENT Cart




Model #
2652LAP
Laparoscopy Cart



Model #
2652ORTH
Ortho Cart

LogiCell® Supply Carts

Triple Cart Shell Configurations



Model #2653LCTS
1 standard cell left
1 wide cell with top split right



LogiCell® Supply Carts

Pre-Configured Triple Carts for OR



Model #
2653VASC
Vascular Cart



Model #
2653ARTH
Arthroscopy Cart



Model #
2653CYSTO
Cysto Cart



Model #
2653CARD
Cardiac Cart



Model #
2653PACK
Drape/Pack Cart



Model #
2653STAP
Staple Cart



Model #
2653ORTH2
Ortho Cart



Model #
2653LAP
Laparoscopy Cart



Model #
2653ORTH1
Ortho Cart

LogiCell[®] Supply Carts

Pre-Configured Catheter Carts



Model # 2652CAT1HG

Cart with 7 slides
Shown with hinged doors ("HG")

Catheter Carts in many styles

- Carts come complete with pictured accessories.
- Dividers included with most baskets, tubs and shelves.
- Customize your own cart using shells from pages 4, 5 and 7 and accessories from pages 24 - 26
- Double carts: 28"D x 40"W x 80"H
- Triple carts: 28"D x 57"W x 80"H



Model #2652CAT4

Supply Cart with 5 shelves,
each with 4 dividers

Model # 2652CAT2

Cart with 7 cathslides and top
shelf with 5 dividers




Model #2652CAT3

Cath box cart with roll-out box
and top shelf with 5 dividers



**Model #
2652CAT5**

Short/long cath cart with
9 slides and top shelf
with 2 dividers



Cath slides on carts have 9 hooks. Each hook holds 15 catheters.
135 catheter capacity per slide. On Shelf slides are adjustable left and right.

**Model #
2652CAT6**
Cath/Supply Cart with
6 cath slides, top shelf
and baskets/tubs



LogiCell® Supply Carts

Pre-Configured Catheter Carts / Catheter Procedure Carts



Model #
2652CAT7HG
5 shelves and 3 cath slides.
Shown w/ hinged doors ("HG")



Model #
2652CAT8
Supply Cart with baskets,
tubs and shelves



Model #
2653CAT1
Cart with 9 slides

Model #2653CAT3
Cart with 3 slides, shelves, tubs,
lockable drawer and baskets



Model #
2653CAT2
Cart with 14 slides,
shelves, tubs and
baskets



Model #
2653CAT4
Cart with 7 slides,
shelves, tubs and
baskets



Model #5140
Catheters, wires, trays, balloons



Model #5120
Catheters, wires, trays, balloons



Model #5130
Catheters, wires, trays, balloons



Model #5110
Catheters only



LogiCell[®] Tambour Door Cabinetry

11000 Series

From Vision to Reality



Our 11000 Series makes the best use of your limited space




- Use crowded corner space with ease
- Tambour Cabinetry does not cramp valuable treatment area
- Tall, floor and wall cabinets in a variety of sizes and configurations



LogiCell® Tambour Door Cabinetry

11000 Series Accessories



Wire baskets, epoxy coated (3", 5" & 8"H)
 • Long and short dividers dividable in 1" increments


Easily slant baskets, tubs and shelves on the cell system



Tubs, ABS plastic (3", 5" & 8"H)
 • Long and short dividers dividable in 4" increments



Drawers, single and double width, 5"H, lockable



Cath Slides (6 or 9 hook)



Shelves with back lip and divider grooves are reversible and can be slanted

LogiCell® Tambour Door Cabinetry

11000 Series Tall Cabinets



Tall Cabinetry

- Frosty white melamine
- Tambour door, light grey with cassette design
- 18" or 26"D with cassette design
- 24", 36" or 42"W
- Most units complete with cell system of plastic side panels (patent applied for) that accommodates a range of accessories (pages 24-26)
- Toe kick recessed base 4"H for cleaning
- Partial sloped top

Model #11112
Shown with accessories



Model #11123LC
Shown with accessories

Model Number	Overall Dimensions (D x W x H)	Interior Dimensions (D x W x H)	Shelf/Cell Configuration	Accessory Size
11112	18" x 27" x 84"	16" x 24" x 78"	Single column with 23 cells	16"
11113	18" x 36" x 84"	16" x 34" x 78"	Single column with 23 cells	16"
11511	26" x 18" x 84"	24" x 16" x 78"	Single column with 23 cells	24"
11523	26" x 36" x 84"	24" x 34" x 78"	Double column with 46 cells	24"
11123LC	18" x 36" x 84"	16" x 34" x 78"	20 Cells at left, 5 x 8"W shelves at right Top shelf 12"H with 5 dividers 4"H	16"
11123RC	18" x 36" x 84"	16" x 34" x 78"	20 Cells at right, 5 x 8"W shelves at left Top shelf 12"H with 5 dividers 4"H	16"
11124LC	18" x 42" x 84"	16" x 40" x 78"	20 Cells at left, 5 x 14"W shelves at right Top shelf 12"H with 7 dividers 4"H	16"
11124RC	18" x 42" x 84"	16" x 40" x 78"	20 Cells at right, 5 x 14"W shelves at left Top shelf 12"H with 7 dividers 4"H	16"

LogiCell® Tambour Door Cabinetry

11000 Series Catheter Cabinets

Let us design



Model #11153CAT2

your catheter



Model #11153CAT1

cabinet for you



Model #11153CAT3

Model Number	Overall Dimensions (D x W x H)	Interior Dimensions (D x W x H)	Number of catheter slides/ Maximum cabinet capacity	Max Package Length	Shelf Configuration
11153CAT1	18" x 36" x 84"	16" x 34" x 78"	7 x 6-hook slides Up to 630 catheter capacity	Up to 58"	Top shelf 34" each with 5 dividers 4"H
11153CAT2	18" x 36" x 84"	16" x 34" x 78"	7 x 6-hook slides Up to 630 catheter capacity	Up to 72"	No top shelf
11154CAT1	18" x 42" x 84"	16" x 40" x 78"	9 x 6-hook slides Up to 810 catheter capacity	Up to 58"	Top shelf 40" each with 7 dividers 4"H
11154CAT2	18" x 42" x 84"	16" x 40" x 78"	9 x 6-hook slides Up to 810 catheter capacity	Up to 72"	No top shelf
11153CAT3	18" x 36" x 84"	16" x 34" x 78"	14 x 6-hook slides Up to 1260 catheter capacity	Up to 31" each space	Two shelves 34" Top Shelf with 4 dividers
11154CAT4	18" x 42" x 84"	16" x 40" x 78"	18 x 6-hook slides Up to 1620 catheter capacity	Up to 31" each space	Two shelves 40" each with 4 dividers 4"H
11513CAT1	26" x 36" x 84"	24" x 34" x 78"	7 x 9-hook slides Up to 945 catheter capacity	Up to 58"	Top shelf 34" with 4 dividers

LogiCell® Tambour Door Cabinetry

11000 Series Wall and Floor Cabinets

Mix storage and workspace with our wall and floor cabinets



Model #11153WC



Model #11113WCR with accessories

Model Number	Overall Dimensions (D x W x H)	Interior Dimensions (D x W x H)	Shelf/Cell Configuration	Access Size (D)
11151WC	18" x 19" x 31"	16" x 17" x 26"	1 melamine shelf	n/a
11153WC	18" x 38" x 31"	16" x 36" x 26"	1 melamine shelf	n/a
11154WC	18" x 42" x 31"	16" x 40" x 26"	1 melamine shelf	n/a
11113WCL	18" x 38" x 31"	16" x 36" x 26"	8 cells at left, 2 x 10" shelves at right	16"
11113WCR	18" x 38" x 31"	16" x 36" x 26"	8 cells at right, 2 x 10" shelves at left	16"
11114WCL	18" x 42" x 31"	16" x 40" x 26"	8 cells at left, 2 x 14" shelves at right	16"
11114WCR	18" x 42" x 31"	16" x 40" x 26"	8 cells at right, 2 x 14" shelves at left	16"

Model #11553FC



Model #11511FC with accessories



Model Number	Overall Dim. (D x W x H)	Interior Dim (D x W x H)	Shelf/Cell Configuration	Access Size (D)	Door Style
11511FC	25" x 19" x 36"	24" x 17" x 25"	8 cells	24"	Tambour
11511FCL	25" x 19" x 36"	24" x 17" x 30"	8 cells	24"	Solid left
11511FCR	25" x 19" x 36"	24" x 17" x 30"	8 cells	24"	Solid right
11553FC	25" x 38" x 36"	24" x 36" x 25"	1 melamine shelf	n/a	Tambour
11553FCS	25" x 38" x 36"	24" x 36" x 30"	1 melamine shelf	n/a	Solid
11523FCS	25" x 38" x 36"	24" x 36" x 30"	16 cells - 8 cells/column	24"	Solid
11523FC	25" x 38" x 36"	24" x 36" x 25"	16 cells - 8 cells/column	24"	Tambour
11554FCS	25" x 42" x 36"	24" x 40" x 30"	1 melamine shelf	n/a	Solid
11536FSK	25" x 38" x 36"	24" x 36" x 30"	Sink unit, no shelf. Sink and faucet not included	n/a	Solid

LogiCell[®] Casework

32000 Laminate/38000 Stainless Steel

LogiCell[®] Casework

A world of size and color possibilities



32000 Series Laminate

32000 Series Laminate available in a wide array of colors, patterns, and finishes.

Custom? No problem
We can fit the most challenging healthcare spaces

38000 Series Stainless Steel


38000 Series Stainless Steel with durable aluminum and glass doors (solid doors available)




LogiCell[®] Casework

Sample Casework Elevations

We work with planners, architects and end users to ensure the optimal use of storage



Sample Cath Lab/IR Elevation



Sample Cath Lab/IR Elevation



Sample OR Elevation



LogiCell® Casework

Tall Cabinets

Contact Customer Service for details on the more than 200 cabinet size options and configurations in casework laminate and stainless steel.



Model #32511LG
shown with accessories 26"D x 19"W



Mount drawers, baskets, tubs and shelves in the modular cell system

- Standard cabinet depths: 18" or 26"
- Standard cabinet widths: 19", 24", 36", 42" or 48"
- Cabinet height: 84"
- Custom widths and heights available
- Cabinets equipped with adjustable melamine shelves and/or cell system
- Cell system of plastic side panels (patent applied for) accommodates a range of accessories (pages 24-26)
- Choose from solid or glass hinged doors
- Toe kick recessed base 4"H for ease of cleaning
- Options include sloped top and pass-through units



Model #32553HG
with 1 fixed middle and 4 adjustable shelves, 26"D x 36"W

LogiCell[®] Casework

Wall and Floor Cabinets

- Floor cabinets: 26" standard depth
- Floor cabinet height: 31"H (36" overall with work surface)
- Standard wall cabinet depth: 13"
- Wall cabinet height: 30"
- Standard cabinet widths: 18", 24", 36", 42" or 48"
- Custom sizes available



Model #32353WCG
wall cabinet with standard shelf 13"D x 36"W

Wall Cabinet Choices

- Solid or glass doors
- Adjustable shelves only



Model #32511FC
shown with accessories



Model #32523FCS
with cell system or #32553FCS with adjustable shelf 26"D x 36"W

Floor Cabinet Choices

- Cell system
(available in 18" and 36" widths)
- Adjustable shelves
- Drawers on full-extension slides
- Solid doors only



Model #32551FCD4

LogiCell[®] Casework

Pre-Configured Catheter Casework



Model #32113CAT1




Model #32113CAT2



Model #32113CAT3

Model Numbers	Overall Dimensions (D x W x H)	Number of catheter slides/ Maximum cabinet capacity	Max Package Length	Shelf Configuration
32113CAT1 38113CAT1	18" x 36" x 84"	7 x 6-hook slides Up to 630 catheter capacity	Up to 58"	Top shelf 34" with 5 dividers 4"H
32113CAT2 38113CAT2	18" x 36" x 84"	7 x 6-hook slides Up to 630 catheter capacity	Up to 72"	No top shelf
32114CAT1 38114CAT1	18" x 42" x 84"	9 x 6-hook slides Up to 810 catheter capacity	Up to 58"	Top shelf 40" with 7 dividers 4"H
32114CAT2 38114CAT2	18" x 42" x 84"	9 x 6-hook slides Up to 810 catheter capacity	Up to 72"	No top shelf
32113CAT3 38113CAT3	18" x 36" x 84"	14 x 6-hook slides Up to 1260 catheter capacity	Up to 31" each space	Two shelves 34" and 5 dividers 4"H
32114CAT3 38114CAT3	18" x 42" x 84"	18 x 6-hook slides Up to 1620 catheter capacity	Up to 31" each space	Two shelves 40" and 7 dividers 4"H
32513CAT1 38513CAT1	26" x 36" x 84"	7 x 9-hook slides Up to 945 catheter capacity	Up to 58"	Top shelf 34" with 5 dividers

LogiCell[®] Casework



38000 Series Stainless Steel Cabinets
with LogiCell cell system for OR, Cath Lab, and IR

38000 Series Stainless Steel

38000NDS Nurse Documentation Station



38000 Series Stainless
Steel Scope Cabinets



LogiCell® Scope Cabinets



Model #32636SCHG-9
with slope top option



H.D. hinges, 5 knuckle hospital grade

Vent portal

10% More Capacity than Competitive Cabinets!



Drip Tray

Antibacterial lower padded wall

LogiCell® Scope Cabinets

- Full body scope racks to protect expensive equipment bodies
- Fully padded lower cabinet back wall, antimicrobial
- Drip Tray, removable & cleanable
- Fully vented cabinets
- Extra capacity over competitors, 9 scopes in 36" vs 8 scopes

Model #32636SCHG-9
with slope top option

LogiCell® Scope Cabinets



38636SCHG9-Custom



11636SC - 9



Double scope depth

MODEL NO.	DIMENSIONS	DESCRIPTION
11624SC-5	18"D x 27"W x 92"H	5 Scope
11636SC-9	18"D x 36"W x 92"H	9 Scope
11636SC-18	18"D x 36"W x 92"H	14 Broncho Scopes, two level
11636SC-18DB	18"D x 36"W x 92"H	18 Full Size Scopes, 1 level dbl. deep
11642SC-11	18"D x 42"W x 92"H	11 Scope
11642SC-22DB	18"D x 42"W x 92"H	22 Full Size Scopes, 1 level dbl. deep
111824ST	18"D x 27"W	Slope Top
111836ST	18"D x 36"W	Slope Top
111842ST	18"D x 42"W	Slope Top

* Full slope top available as accessory when specified, 96"H

11000 Series

- Tambour door
- Frosty white melamine
- Toe kick base, 4" recess
- Drip Tray, removable & cleanable
- Partial slope top *

32000 Series

- Glass or solid doors
- Frosty white melamine
- Toe kick base, 4" recess
- Drip Tray, removable & cleanable
- Flat top *


38000 Series

- Stainless cabinets
- Same sequence of sizes & options as 32000 above

MODEL NO.	DIMENSIONS	DESCRIPTION
32624SCLG-5	18"D x 27"W x 92"H	5 Scope, glass door, left hinge
32624SCRG-5	18"D x 27"W x 92"H	5 Scope, glass door, right hinge
32624SCLS-5	18"D x 27"W x 92"H	5 Scope, solid door, left hinge
32624SCRS-5	18"D x 27"W x 92"H	5 Scope, solid door, right hinge
32636SCHG-9	18"D x 36"W x 92"H	9 Scope, glass door
32636SCHS-9	18"D x 36"W x 92"H	9 Scope, Solid door
32636SCHG-18	18"D x 36"W x 92"H	18 Broncho Scopes, TwoLevel, Glass Door
32636SCHS-18	18"D x 36"W x 92"H	18 Broncho Scopes, TwoLevel, Solid Door
32636SCHG-18DB	18"D x 36"W x 92"H	18 Full Size Scopes, OneLevel, Dbl Deep, Glass Door
32636SCHS-18DB	18"D x 36"W x 92"H	18 Full Size Scopes, OneLevel, Dbl Deep, Solid Door
32642SCHG-11	18"D x 42"W x 92"H	11 Scope, Glass Door
32642SCHS-11	18"D x 42"W x 92"H	11 Scope, Solid Door
32000ST		Slope Top, Sold by the linear foot

LogiCell® Accessories

Baskets and Shelves




Easily slant baskets, tubs and shelves on the side panels




Model #LH-3 Label Tag
in white, yellow, green, blue & red

Wire Basket Advantages


- 3" and 5" baskets slant easily without added parts
- High density: Dividable in 1" increments
- Integrated safety stop
- Full access, 80% extension
- Epoxy coated steel fabrication
- 24" x 16"



Model #3B Basket 3"H
Occupies 1 cell
Long Divider #3BDL
Short Divider #3BDS



Model #5B Basket 5"H
Occupies 2 cells
Long Divider #5BDL
Short Divider #5BDS



Model #8B Basket 8"H
Occupies 3 cells
Long Divider #8BDL
Short Divider #8BDS



Shelves accommodate cath slides and dividers



Shelf slanted with lip in front

Shelf Advantages

- Slants easily
- Lip can be used in back or front
- Safety stop prevents shelf from falling out
- Numerous divider grooves
- Pre-engineered to hold catheter slide ("SHC")
- Epoxy coated steel fabrication
- Occupies 2 cells with divider

Model #1624SHC 16"D shelf for single column
Model #1636SHC 16"D shelf for double column
Model #164DIV Divider 4"H

Model #2416SHC 24"D shelf for single column
Model #2436SHC 24"D shelf for double column
Model #244DIV Divider 4"H

Economy Shelves

- Slants easily
- Can accommodate dividers but not cath slides
- 35 lb. capacity. Occupies 2 cells with dividers

Model #1624SHE 16"D shelf for single column
Model #2416SHE 24"D shelf for single column



Model #2416SHC with divider



Model #244DIV divider

LogiCell® Accessories

Tubs and Drawers

Plastic Tub Advantages

- 3" and 5" tubs slant easily without added parts
- High density: Dividable in 4" increments
- Integrated safety stop
- Full access, 80% extension
- White, ABS plastic
- 24" x 16"

Fully extended basket and tubs



Model #3T Tub 3"H

Occupies 1 cell
Long Divider #3TDL
Short Divider #3TDS



Model #5T Tub 5"H

Occupies 2 cells
Long Divider #5TDL
Short Divider #5TDS



Model #8T Tub 8"H

Occupies 3 cells
Long Divider #8TDL
Short Divider #8TDS

Drawer Advantages

- Frosty white laminate
- Heavy-duty full extension slides
- Lockable
- Pre-mounted inside cabinet at requested position
- Easily change location
- Top area functions as shelf
- Occupies 2 cells, 3.85" interior height

Single column drawer 18"D



Double column drawer 18"D



Model #1624LDR 16"D drawer for single column
Model #1636LDR 16"D drawer for double column


Model #2416LDR 24"D drawer for single column
Model #2436LDR 24"D drawer for double column

Single and double column drawers 24"D



Accessories

Catheter Accessories



LogiCell® Catheter Slide for 16"D Shelf

- Model #CAT16
- Epoxy powder coated steel
- 6 hooks
- Each hook can hold up to 15 packages
- Easy pull mechanism
- Includes label holder at front

LogiCell® Catheter Slide for 24"D Shelf

- Model #CAT24
- Epoxy powder coated steel
- 9 hooks
- Each hook can hold up to 15 packages
- Easy pull mechanism
- Includes label holder at front

Wire Shelf Stainless Steel Catheter Slides

- Mount beneath any wire shelf
- Card holder included
- Each hook holds up to 15 packages
- **Model #USC24** 24"D slides with 7 hooks
- **Model #USC18** 18"D slides with 5 hooks



Model #CHC Catheter Hanging Clips

- Adhesive bottom attaches to non-punched or punctured packages
- Hanger hook at top
- Available in packages of 180

Model #CSM Catheter Stock-Manager Clips

- Attaches to the end of cath slide hooks
- Available in packages of 250

Model #CHB Catheter Hanging Bar

- Available up to 36"L
- Capacity of approx. 7 hanger clips

Wire Products

Functional, Economical

and Adaptable...Wire works well!

Wire Catheter and Supply Carts



Model #836CH

Model Number	Overall Dimensions (D x W x H)
836CH	18" x 36" x 63"
856CH	18" x 48" x 63"
866CH	18" x 60" x 63"
456CH	24" x 48" x 63"
466CH	24" x 60" x 63"
476CH	24" x 72" x 63"

Wire Shelving Units

- Nickel chrome plated steel
- Complete with 4 shelves adjustable on 1" increments
- Add prefix "MB" to model # to add 5" poly swivel casters two with brakes (also adds 6" to height)
- Many other accessories available such as dividers, ledges, label holders, end panels, bins, totes
- Dozens of other standard sizes and configurations available

Model #CATCOMB2CH



Model # CATCOMB3CH

Model # HC2436E-CH



HC Style Catheter Cart with Stainless Slides

Model #HC2436E-CH
 (24" x 36") shown economical catheter storage
 24" (7-hook) or 18" (5-hook) depths
 widths from 24" to 72"

Combination Catheter Cart Shelves and hooks for a range of supply storage

4-level step down catheter rack with 40 hooks, 2 shelves 60"W, 3 shelves 24"W

Model #CATCOMB2CH Open unit 18"D x 60"W x 79"H overall

Model #CATCOMB3CH Same as above but with ends and back and 7 bins

Model #CRK40 also available - Unit with 40 catheter hooks only, 18"D x 36"W x 64"H overall



Stem Caster Cart
 Model #MB456CH



Stem Caster Basket Cart

- 18"D x 36"W x 78"H
- 6 each baskets
- 24 each dividers
- 2 - 18" x 36" shelves
- 5" casters

Model #AF090CH
 with DV2010CH

Wire Products

IV Carts



26517IV



500IV-III

500IV-III

- Open wire cart
- All welded frame
- 4 - 8", 1 - 5" basket
- Top basket shelf
- Solid base
- 6" casters

#27517IV

- Closed, secure cart
- Key Lock
- 4 - 8", 2 - 5" baskets
- Ergonomic handles
- 6" casters



Model #26517
IV Cart



Full Extend Baskets



Model #500IV-III
IV Cart




Full Extend Baskets

Wire Products

Par Wall[®] System "One Product, One Lane"

- Supplies are organized and accessible. No more "digging"
- Choose from wall-mount louvered panels or free-standing epoxy-coated frame
- Regular or fine mesh wire baskets in 12" or 18" depths and 18" or 36" widths
- Dividers go front to back. Smaller dividers for fine mesh create sub-divided spaces for small supplies
- Contact your space solutions expert for a consultation




Divide front to back or side to side



Both systems use the LH-3 label holder



Par Wall[®] Frame with regular and fine mesh baskets



Wall mount style with regular baskets

ParStor-III[®] Ultra High Density Storage

- Increases your storage capacity 30-45%
- Opens up 30-50% more floor space
- Open wire minimizes dust accumulation
- All baskets/drawers are on premium full-extension slides
- Standard module includes 5 - 5"H baskets, 1 - 5"H drawer, 1 - 3"H baskets and 1 - 8"H basket. Top wire shelf, bottom poly insert and 3" swivel casters, two with brakes.
- Components easily configured to your needs
- All welded Par Carts are also available, like the #500IV-III on the facing page



All frames are numbered



Long and short dividers



Safety stop on slanted baskets

Samples, In Use

Closed System



LogiCell®

Open System



ParStor-III®

Interchangeable Accessories



Cabinets in custom finish laminate to match any decor



Custom design storage in really odd corners!
Use all your space, don't let the space use you.



Open Par Wall® Secure Closed LogiCell®

Samples, In Use



Cath Lab Core



Cath Lab Core, closed tambour doors



11000 Series, tambour door cabinets in a wide range of sizes and a full range of accessory options

Samples, In Use



Sterile Clean Storage Room



Full Extend Accessories



Glass Doors

Stainless Steel Scope Cabinet



Stainless Steel Scope Racks

Terms & Conditions

Acceptance

All orders subject to acceptance of Terms and Conditions of sale, regardless of the form in which submitted. All orders are quoted F.O.B. Galesburg, Michigan unless stated otherwise. Custom manufactured products may be subject to a deposit at the time of order.

Order Confirmation

To avoid duplication of shipment, the word "CONFIRMATION" or "CONFIRMING" should appear on all documents which are sent to confirm orders placed by telephone, fax or any other way.

Terms

Minimum order \$50.00. Terms are Net 30 days (subject to credit approval). An annual service charge of 24% will be assessed on past due accounts.

Shipments

Shipments will be made F.O.B. shipping point and the buyer assumes responsibility for all damages in transit. Such damages must be declared at time of receipt, directly on the shipping documents, and claims must be filed immediately with the transport company. If after receipt of shipment, concealed damage is discovered, immediately notify the delivering carrier requesting an inspection. The inspection report should be used to support your claim against the carrier. Any claim for damage or shortage must be filed by you with the delivering carrier. Unless a preferred routing is specified at the time of order, shipment will be made the best way, transportation prepaid and added to invoice.

Returned Goods Policy

1. All returns must be authorized.
2. Each request must include the following:
 - A. Customer name and address.
 - B. Invoice number
 - C. Invoice date
 - D. Customer P.O. number
 - E. Quantity, catalog number and description of item
 - F. Reason for return
3. Upon receipt of your request, a return authorization number will be issued and this must be clearly referenced with your return.
4. The return shipment must be coordinated with a customer service representative, who will advise the proper mode of transport for the return.
5. Returned goods may be subject to a minimum 20% handling charge based on state, condition of goods and reason for return.

Returned Goods Policy

6. Returns will NOT be authorized for the following:
 - A. Items purchased on a "special order" basis (unless original vendor will accept the return).
 - B. Items that are custom manufactured.
 - C. Partial units of sale.
 - D. Items without original packaging materials.
7. **RETURNS MUST BE PREPAID AND WITHIN 10 DAYS OF ORIGINAL SHIPMENT.**

Specifications

All dimensions and shipping weights are approximate. LogiQuip reserves the right to change specifications or substitute products of equal specifications without notice. The company also reserves the right to cease manufacturing or distribution of any product without notice. Manufacturer's specifications and prices are subject to change without notice.

Freight Classification

Freight Class is product specific and may vary from CL70 to CL250.

Installation

Installation of cabinets: Must be in compliance with all state and local code regulations and should be performed by experienced and qualified tradespeople. For best long term results and proper functioning of the cabinet doors and latches, the cabinets must be installed plumb and squarely true with adjoining cabinets. Compliance with all state and local regulations relative to anchoring cabinets should be confirmed by the appropriate hospital staff. Improper installation, or failure to comply with state and local regulations, may result in failure of the product and denial of potential warranty claims.

Warranty Statement

LogiQuip equipment is warranted to be free from defective material and workmanship under normal use and service for a period of one year from date of shipment to customer. During this period, LogiQuip will repair or replace (at its option) without charge, F.O.B. plant of manufacture, the equipment or any part thereof which factory inspection discloses to be defective. Shipments of defective material for repair must be sent prepaid.



LOGICELL®

By **+** LOGIQUIP

Phone: 800.665.3760 | www.LogiQuip.com
Fax: 888.965.3984 | info@logiquip.com