



umano**medical**

## LONG TERM CARE BED WITH BED EXIT DETECTION SYSTEM



10" low height with built-in bed exit detection system to minimize falls

Seamless design to maximize infection control

Four casters with central brake to maximize staff efficiency and patient safety

## CLEANABILITY • SCALABILITY • LOW HEIGHT SAFETY



Clearance for lift and overbed table compatibility, even in the lowest height



Four casters with central brake to minimize falls



Built-in bed exit system for improved detection reliability  
Available with or without display screen



Optional footboard design to fit your environment

## LONG TERM CARE BED WITH BED EXIT DETECTION SYSTEM

### BED FEATURES

- Built-in bed exit alarm system
  - 2 levels sensitivity
- 4 DC Motor
- Battery back-up
- 4 casters with central brake
- Flat deck for ease of cleaning
- Open architecture
- One step Trendelenburg
- Patient control lockout
- Boostless™ backrest system
- Compliant with standard IEC 60601-2-52

### OPTIONAL FEATURES

- Built-in bed exit alarm system
  - 3 levels of sensitivity (with Smart screen)
- Remote bed exit alarm connectivity
- Auto-contour
- Backrest angle indicator
- One step vascular foot position
- Programmable lowest height + 2" / 5 cm
- 80" / 203.2 cm patient sleep surface
- 84" / 213.4 cm patient sleep surface
- 5" / 12.7 cm single or dual casters
- 120-volt/5A insulated outlets with USB mobile device power outlets (not available for 220 V)
- CPR release handle
- Dual nightlights
- Screen on footboard control
  - Bed position indicator
  - Bed exit events log
  - Inform Message center
  - Maintenance menu
  - Password management
- 2 functions patient pendant control
- 5 functions patient pendant control
- Plastic head & footboard
  - Transport handles
  - Integrated pump holder
  - Optional display screen
- Wooden head & footboard
  - Choice of customized colors
- Patient mobilization helper
- Foot end back-up control on the frame
- Patient control connector (Permanent or removable)

### Accessories\*

- Patient helper system
- Upright oxygen bottle holder
- IV poles (Permanent & Removable)

\* To see all accessories available, visit our website

### Specifications for model FL36

Overall length	88" / 224 cm
Overall width	
Siderails down	39.25" / 100 cm
Siderails up	40.5" / 102.9 cm
Weight capacity	500 lb / 227 kg
Height range (to litter top)	
Low (with 5" / 12.7 cm dual casters)	10" / 25.4 cm
High (with 5" / 12.7 cm dual casters)	30" / 76.2 cm
Low (with 5" / 12.7 cm single casters)	11.5" / 29.2 cm
High (with 5" / 12.7 cm single casters)	31.5" / 80 cm
Mattress deck positioning	
Backrest	0° to 60°
Knee Gatch	0° to 30°
Trend. /Reverse Trend.	± 16°
Boostless™	4" / 10.2 cm
Patient surface	35" x 80" or 84" 88.9 x 203.2 or 213.4 cm
Casters diameter	5" / 12.7 cm
Bed weight	336 to 369 lb / 152 to 167 kg

### Electronics

Compliant with standard	IEC 60601-1 3 <sup>rd</sup> edition
Volts	120 V AC / 220-240 V AC
Ampere rating	4,25 A (9,25 A with insulated outlets) / 2,15 A
Current leakage	< 300 micro amps
Frequency	60 Hz / 50 Hz
Hospital grade plug	Yes
Battery type	Sealed rechargeable lead acid 1.2 Ah, 2 * 12 V

*Umano Medical* reserves the right to change specifications without notice.  
Contact your local representative to know which configuration is available in your region.



### Umano Medical

230, boul. Nilus-Leclerc, L'Islet (Québec) G0R 2C0 CANADA

Toll free number : 1 844 409-4030

information@umanomedical.com

www.umanomedical.com