



ook«snow»

LOW HEIGHT • CLEAN • EASY TO USE • FULLY FEATURED

ATTENTION TO NEEDS

There are many elements to medical care, however the majority of them evolve around the bed.

A dignified medical environment that offers appealing beds goes a long way towards easing their transition from home to hospital.

Umano Medical beds offer many options to healthcare professionals and can be adapted to meet the requirements of a demanding hospital environment. Our ergonomic designs give patients more autonomy and ensure their safety.



CARE
GIVER



"No matter how busy I am, I give my patients personal attention and I involve them in all aspects of their care."

www.umanomedical.com

INTUITIVE

EASY TO USE

Umano Medical beds are easy to operate with minimal training. Patients can easily access the control panel to position the bed, reducing the frequency of calls for assistance and allowing them to feel more in control of their surroundings.

Our advanced technological features have been designed to deal with any medical emergency. Commands are accessible through easy to recognize visual prompts, providing healthcare professionals with enhanced quick-response capacity and allowing them to dedicate more time attending to the patient's needs.



MOBILITY



One-hand
siderail release



Central brake
pedals



One-hand
emergency
CPR release

LOW HEIGHT HOSPITAL BED BUILT WITHOUT COMPROMISE

MEDICAL & SURGERY BEDS WITH SCALE AND INTEGRATED BED EXIT SYSTEM


The **ook«snow»** Med Surg Bed offers best in class low height without compromising important features such as lift and overbed table compatibility, one-step Trendelenburg and braking.

Designed to maximize infection control and reduce falls, **ook«snow»** improves patient's well-being thanks to its fully featured design, its built-in safety characteristics and its ease of use.

ook«snow» reduces bed inventory and purchasing costs by combining both low bed and hospital bed characteristics.



Different graphics are available to suit your facility.



Clearance for lift and overbed table compatibility (5" under bed clearance)



Open architecture for ease of cleaning and simple installation of mattress



Boostless™ backrest system for improved patient comfort during bed activation



Optional integrated bed extender for more flexibility

FULLY FEATURED



USB POWER OUTLETS

Conveniently located USB ports allowing patients to maintain their mobile devices to a full charge.



FALL RISK MANAGEMENT

Keep your patients safe with the right combination of low height, Smart Bed Exit technology and ergonomic siderail design.



10" LOW HEIGHT

A low bed that reduces the adverse effect of a fall without compromising other important characteristics.



PATIENT CONTROL LOCK

Allows only safe bed movements.



FOOT BOARD CONTROL PANEL

Advanced caregiver controls with Smart Screen to communicate important safety messages or access a log of patient weights and bed exit events.



USER-FRIENDLY CLEANING ARCHITECTURE

Seamless smooth surfaces of an open bed architecture eases infection control.



BUILT-IN EFFICIENT SCALE

Smart scale providing accurate weight at any bed position.



3-STEP BED EXTENDER

Built-in bed extender increases patient comfort and also caregiver efficiency.



Central brake



Bed exit side view lights



Back-up control on the frame



USB power outlets



3-step bed extender - 80" - 84" - 90"

SIMPLIFY CLEANING

We know that properly cleaned and sanitized beds are essential to infection and disease control in hospitals. The responsibility to ensure the highest expectations rests by Environmental Service Workers (ESWs). At *Umano Medical*, we believe in sharing that same responsibility by designing our beds to make their job easier with the following features:

- no special tools needed for cleaning
- seamless smooth surfaces and finishes
- siderail bottoms that never touch the floor
- flat mattress deck
- easy access to all parts
- hidden cables and wires
- removable or replaceable deck covers



USER-FRIENDLY CLEANING ARCHITECTURE

Seamless smooth surfaces of an open bed architecture eases infection control.



FALL RISK PREVENTION

INCREASED SAFETY, FOR YOU AND YOUR PATIENTS

Promoting and encouraging patient autonomy and mobility is an important element in providing hospital care. Whether they are in a short-term or a long-term care facility, patients have an understandable fear of falling. For the healthcare provider, the risk of serious injury caused by forceful exertions and non-ergonomic positions is also a concern.

Umano Medical beds are designed and built with that awareness.

Our beds offer affordable options and features for every scenarios:

- **Low height** allows smaller, frail and post-surgery patients to benefit from the safety of a 10" low height bed.
- **A local or remote alarm** can detect patient movements and can be adjusted according to the level of risk.
- **Built-in bed exit and nightlights** allow fall detection reliability while encouraging patient autonomy. Under bed lights also provide patient and staff with a clear view of the floor and around the bed.
- **One-touch adjustment features** make safe bed positioning an easy task.



HEALTHCARE
PROVIDER



" I am always looking out for the possibility that my patients may fall out of bed. For me, a medical bed that allows me to have more control over their safety goes a long way toward giving me peace of mind "

www.umanomedical.com



Toll free number : 1 844 409-4030

www.umanomedical.com