



bassinets





shown with optional bassinet and pad

P-1110-SS bassinet stand

frame size

18 $\frac{1}{4}$ "W x 31"L x 37 $\frac{3}{4}$ "H

Height with bassinet basket - 41"H

P-1110

- Same as P-1110-SS except frame is powder coated steel

features

- Design allows close approach from both sides
- Basket can be raised at either end
- 3" double ball bearing swivel casters
- Easy to clean
- Stainless steel construction
- *Lower shelf available as an option*



shown with optional bassinet and pad

P-1110-A-SS bassinet stand with drawer

frame size

18 $\frac{1}{4}$ "W x 31"L x 37 $\frac{3}{4}$ "H

Height with bassinet basket: 41"H

P-1110-A

- Same as P-1110-A-SS except frame is powder coated steel

features

- End mounted drawer measures: 16 $\frac{1}{2}$ "W x 13 $\frac{1}{4}$ "L x 7 $\frac{3}{8}$ "H
- Design allows close approach from both sides
- Basket can be raised at either end
- 3" double ball bearing swivel casters
- Easy to clean
- Stainless steel construction
- *Lower shelf available as an option*



shown with optional bassinet and pad

P-1110-B-SS bassinet stand with drawer

frame size

18 $\frac{1}{4}$ "W x 31"L x 37 $\frac{3}{4}$ "H

- Height with bassinet basket: 41"H

features

- Side mounted drawer measures: 16 $\frac{1}{2}$ "W x 21 $\frac{1}{4}$ "L x 3 $\frac{1}{8}$ "H
- Design allows close approach from both sides
- Basket can be raised at either end
- 3" double ball bearing swivel casters
- Easy to clean
- Stainless steel construction
- *Lower shelf available as an option*



shown with optional bassinet and pad

P-1114-SS bassinet stand with tote box drawers

dimensions

32"W x 20"L x 33"H

Height with bassinet basket: 41"H

features

- Design consistent with rooming in techniques
- Two 10" deep tote boxes accessible from either side
- 4" double ball bearing swivel casters
- Easy to clean
- Stainless steel construction



shown with optional bassinet and pad

P-1115-A-SS bassinet cabinet

dimensions

19 $\frac{3}{4}$ "W x 30 $\frac{7}{8}$ "L x 37 $\frac{1}{4}$ "H

* *Bassinet not included.*
Height with bassinet 44 $\frac{1}{8}$ "

Cabinet dimensions

17"W x 29 $\frac{1}{2}$ "L x 24 $\frac{7}{8}$ "H

Drawer frame with removable pan

11 $\frac{1}{4}$ "W x 16 $\frac{1}{2}$ "L x 4"H

features

- Bassinet carriage slides back to allow working area
- Friction lock secures the bassinet in both open and closed positions
- Removable internal shelf, positioned in middle left compartment
- 4" double ball bearing swivel casters
- Easy to clean
- Stainless steel construction

P-1115-SS

Same as P-1115-A-SS except without doors

options

Bassinet Basket (Part No. 53) Basket measures 16"W x 30 $\frac{1}{4}$ "L x 9 $\frac{5}{8}$ "H

Bassinet Pad (Part No.54) Pad measures 12 $\frac{1}{2}$ "W x 26 $\frac{1}{2}$ "L x 1" Thick

