

## Urology Drain Bags



### Square Drain Bag

**BRA5417**

with plastic tab injection molded screen and funnel, Case/20  
Fits all square frames  
Sterile, Individually pouched



### Square Drain Bag

**BRALG2500**

with metal tabs and mesh screen, Case/20  
Fits all square frames  
Sterile, Individually pouched



### Hoop Drain Bag

**BRALG2800**

Fits GE/OEC Uroview 2600/2800 Tables, Case/20  
Heavy gauge vinyl strip holds bag snugly to the foot of the table  
Mesh screen, flexible hose extends 6 feet  
Under buttocks flap  
Sterile, Individually pouched



### Urology Drain Bag

**BRA5442**

Fits Siemens Uroskop System  
Heavy gauge vinyl strip holds bag snugly to the foot of the table  
Mesh screen, flexible hose extends 6 feet  
Under buttocks flap  
Sterile, Individually pouched

**NOTE:** All drain bags are also available as non-sterile.

When ordering include N/S at the end of the product identifier to indicate non-sterile product is required. E.G. BRA5442 N/S