



## Med Cart Canada Inc.

5466B Timberlea Blvd, Mississauga

Ontario, Canada, L4W 2T7

Phone: 1-855-212-9696

Local: 905-212-9696

Fax: 905-212-9692

sales@medcartcanada.com

www.medcartcanada.com

### Smart Flex Cart #SFC40

Dual Column Cart, designed for all bin types

Engineered eMAR ready



#### Construction:

- Painted lightweight aluminum cabinet and drawers

#### Standard Accessories:

- Internal narcotics box with Illinois Lock
- BEST® brand lock on cabinet
- Painted pull-out side shelf with recessed work area
- Ball-bearing, full extension drawer slides
- Plastic replaceable top with durable work surface
- Ergonomic top rail/push handle
- Easy clean protective pontoon bumpers
- 5" ball-bearing premium rubber casters, all full swivel, 1 brake and 1 directional
- Hard-baked powder coat finish (Customer chooses color for cabinet)
- Standard accessories also include aluminum waste container and three section cup and straw holder shown below



#### Thirteen available color choices:

Beige, Burgundy, Light Blue,  
Dark Blue, Light Gray, Mauve,  
Red, Teal, Yellow,  
Hammerstone Grey

Hammerstone Blue, Hammerstone Green and Sand

Colors may be viewed at

[www.medcartcanada.com](http://www.medcartcanada.com)

#### Drawer Configuration:

Right column:

- 4 ea – 4" drawers
- 1 ea – 11" multipurpose storage drawer

Left column:

- 4 ea – 4" drawers
- 1 ea – 11" multipurpose storage drawer with locking narcotics box

#### Assembled Cart Dimensions:

- 42" H x 36" W x 24" D (includes top rail)  
(106.6 cm x 91.4 cm x 61 cm)
- Work surface: 38" H (96.5 cm)