



umanomedical

OVER BED TABLE
FTG® SERIES



**EASING
INFECTION
CONTROL**

Will meet the highest expectations
for infection control

A low profile base for an easy
access under the bed

Extendable acrylic top for increased
reach and greater surface
(optional)

Fast and easy adjustment

**LOW HEIGHT • EASING INFECTION CONTROL
• COMPATIBLE WITH MOST CHAIRS**



Resistant acrylic top designed for ease
of cleaning and disinfection



31" / 78.7 cm wide base
for compatibility
with most chairs



Best in class low height



Umanomedical

**OVER BED TABLE
FTG® SERIES**

STANDARD FEATURES

- Twin Caster 1.97" / 5 cm
- Load Capacity 100 lb / 45 kg
- Self-locking Hi-lo system
- Sturdy construction
- Epoxy polyester powder coating


OPTIONS

CHOICE OF LAMINATES COLOR*




* Possibility of customized colors, contact your local representative

CHOICE OF ACRYLIC TOP COLOR



CHOICE OF BASE COLOR



Specifications

TOP	LAMINATED	ACRYLIC	EXTENDABLE ACRYLIC
Height			
Maximum	40.83" 103.8 cm	40.6" 103.1 cm	40.6" 103.1 cm
Minimum	25.83" 65.6 cm	25.6" 65 cm	25.6" 65 cm
Load Capacity			
	100 lb 45 kg	100 lb 45 kg	100 lb 45 kg
Weight			
	38 lb / 17.3 kg	30 lb / 13.6 kg	30 lb / 13.6 kg
Top size			
	16" x 36" 40.6 cm x 91 cm	17.5" x 35" 44.5 cm x 89 cm	17.5" x 35" 44.5 cm x 89 cm
Top Extension			
	NA	NA	8" / 20.3 cm
Base			
	Painted	Painted	Painted
Height	2.875" / 7.3 cm	2.875" / 7.3 cm	2.875" / 7.3 cm
Width	31" / 78.7 cm	31" / 78.7 cm	31" / 78.7 cm

Umano Medical reserves the right to change specifications without notice.



Umano Medical inc.

230, boul. Nilus-Leclerc, L'Islet (Québec) G0R 2C0 CANADA
 Toll free number : 1 844 409-4030
 information@umanomedical.com
 www.umanomedical.com