

Welch Allyn Scale-Tronix® Pediatric Scales



Your littlest patients require specialized weighing solutions. The Scale-Tronix line of pediatric scales is used for measurement of neonates through toddlers, helping you better care for these sometimes difficult-to-weigh patients.

With increments of one gram, Welch Allyn Scale-Tronix Pediatric Scales even help you measure breastfeeding intake, urine output and a number of other specialty weighing needs.

Providing durability and resolution in a big way for your littlest patients—it's what has made **Scale-Tronix a trusted brand for healthcare environments.**



Durable. Long Lasting. Designed for Demanding Healthcare Environments.



4802D Pediatric/Infant Scale

- Designed for daily and serial weighing of neonates, infants and toddlers
- Displays actual weight, rather than averages, to support proper care and treatment
- With increments of one gram, you can measure breastfeeding intake and urine output

4302 Diaper and Organ Scale

- With increments of one gram, you can measure urine output
- Calibration in grams provides one-to-one relation to cc's of urine for easy, rapid calculations (Pounds and ounces may also be selected.)

Find the Right Scale for You

	4802D Pediatric/Infant Scale	4302 Diaper and Organ Scale
Weighing Capacity	44 lb / 20 kg	Standard Capacity: 3,000 gm Extended Capacity: 5,000 gm
Platform	Overall 27.5" w x 16" d x 7" h stainless steel weighing platform. Removable acrylic cradle: 25" w x 14" d. Length measure in centimeters and inches.	8.5" x 8.5"
Resolution	0.2 oz/5 gm. Diaper/preemie weighing (optional): 1 gm	1 gm
Readout	Digital LED display in pounds and kilograms, or kilograms only	Digital LED in grams, pounds or ounces. Automatic shutoff
Power Source	Battery powered (6 D-sized alkaline batteries). Batteries and power adapter included.	6 C-sized alkaline batteries – 10,000 weighings. Power adapter available
Optional Accessories	Printer (external), mobile weighing cart with dispenser for scale liners and storage shelves, table-top dispenser for scale liners, scale liners (500/carton), battery pack	Power adapter available
Automatic Zero	Automatically returns to zero, ready for next patient	N/A
Weight Recall	Recalls last weight	N/A
Reweigh	Recomputes patient's weight while patient is on scale	N/A
Computer Capability	RS-232 output (optional)	N/A

Welch Allyn, Inc.
4341 State Street Road
Skaneateles Falls, NY 13153 USA
(p) 800.535.6663 (f) 315.685.3361