

Wall-Mounted Tubular Rack

Chrome Lead Apron and Glove Holder, Holds one apron and one pair of gloves. No assembly required. Installation hardware is not included.



Code	Details	Dimensions	Weight
ART	Wall-Mounted Tubular Rack	27.75"W x 17"H x 5 ¹ / ₄ "D	8.5 lbs

Wall-Mounted Compact Rack

Ideal wall mount rack for storing multiple pairs of gloves and aprons. Also has convenient bar to hold thyroid collars. Sturdy steel construction. No assembly required. Installation hardware is not included.



Code	Details	Dimensions	Weight
ARC	Wall-Mounted Compact Rack	19" W x 20" H x 8"D	10lbs

Wall Mount Full Swing Rack

Excellent storage capacity. This heavy duty metal rack holds up to five aprons on 5 pivoting apron hangers. Each hanger can swing a full 180 degrees to either left or right. Maximum suggested weight capacity of 110 pounds. Some assembly required. Installation hardware is not included.



Arms perpendicular to the wall: 22"

Code	Details	Dimensions (Plate)	Weight
AR6FS	Wall-Mounted Left and Right Folding Apron Rack	17" W x 4 ⁵ / ₈ " H	15lbs

Wall Mounted Towel Bar Rack 1 & 2

When mounted securely to a wall, can hold one apron (but will accommodate more). 1.75" diameter bar will prevent apron creasing and possible lead damage. Glove holders are mounted onto back panel and located side by side. Available with glove holder for 1 pair or 2 pairs of gloves. All steel heavy duty construction. No assembly required. Installation Hardware is not included.

Code	Details	Dimensions	Weight
ARTB1	Wall Mounted Towel Bar Rack	28" W x 7" H x 6 ¹ / ₂ " D With glove holder 25" H	13lbs
ARTB2	Wall Mounted Towel Bar Rack	28" W x 7" H x 6 ¹ / ₂ " D With glove holder 25" H	14lbs



Wall-Mounted Peg Rack

Peg racks are a quick, easy, durable storage solution. They are constructed of steel with an electro-magnetic powder coating. When mounted securely to a wall, peg racks can hold anywhere from 1-9 aprons depending on size. All pegs measure 4.5" L and are slanted to prevent apron from slipping off peg. No assembly required. Installation Hardware is not included.

Code	Details	Dimensions	Weight
ARP1	1 Apron / 2 peg	6"W x 3"H x 4.5"D	4.2lbs
ARP2	2 Aprons / 4 peg	12"W x 3"H x 4.5"D	6lbs
ARP4	4 Aprons / 8 peg	19"W x 3"H x 4.5"D	12lbs
ARP7	7 Aprons / 14 peg	36"W x 3"H x 4.5"D	13.6lbs
ARP9	9 Aprons / 18 peg	50"W x 3"H x 4.5"D	18lbs





Five Arm Lead Apron Wall Rack

Holds up to five aprons on five arms. It features five chromium-plated pivoting hangers which can flatten out nearly 180 degrees for space saving and storage behind doors, equipment, and more. Arms can swing either left or right. Weight capacity 600lbs. No assembly required. Installation hardware is not included.

Code	Details	Dimensions	Weight
AR6L	Wall-Mounted Left to Right Apron Rack	Arms flush against the wall 28"W x 9" H x 4 ¹ / ₈ " D	15lbs
AR6R	Wall-Mounted Right to Left Apron Rack	Arms flush against the wall 28"W x 9" H x 4 ¹ / ₈ " D	15lbs



Mobile Maxi-Rack with AR6L and AR6R

The mobile Maxi-Rack is the solution when wall mounting AR6 Racks are not possible or mobility is required. Supports 1 left folding and 1 right folding rack. Holds 10 aprons on pivoting arms. Weight capacity 600lbs. Equipped with 4 swivel castors (2 with brakes). Some assembly required.

Code	Details	Dimensions	Weight
ARMX	Maxi-Rack complete with 1 AR6L and 1 AR6R	54"W x 54"H x 20"D	125lbs



Vest & Skirt Mobile Rack

Heavy-duty metal apron rack with locking casters. Holds 10 vest and skirt aprons.

Code	Details	Dimensions	Weight
ARVS	Vest & Skirt Mobile Rack	39"W x 55"H x 26"D	120lbs



Mobile Pedestal Apron Rack

Pedestal design is counter weighted for the utmost stability. Holds up to 10 aprons on pivoting arms that fold flat. Will safely hold the heaviest lead.

Code	Details	Dimensions	Weight
ARMX-HB	Mobile Pedestal Apron Rack	20"W x 49 ^{1/2} "H x 24"D	60lbs

Mobile Apron Rack

Heavy-duty metal apron rack with 4 swivel casters (2 casters are equipped with brakes for locking). Hangers not included.

Small: holds 5-8 aprons or a total weight of 150 lbs.

Large: holds 10-13 aprons or a total weight of 300 lbs.

Code	Details	Dimensions	Weight
AR5	Small Mobile Apron Rack	23 ^{1/2} " L" 54"H x 22 ^{1/2} "W	28lbs / 32lbs w hangers
AR10	Large Mobile Apron Rack	35 ^{1/2} " L" 54"H x 22 ^{1/2} "W	33lbs / 41lbs w hangers

KITS


Code	Details
AR5KIT-O	Holds 5 to 8 Aprons includes 5 hangers - OPEN HANGERS
AR5KIT-C	Holds 5 to 8 Aprons includes 5 hangers - CLOSED HANGERS
AR10KIT-O	Holds 10 to 13 Aprons includes 10 hangers - OPEN HANGERS
AR10KIT-C	Holds 10 to 13 Aprons includes 10 hangers - CLOSED HANGERS



Apron Hanger

A chrome plated steel hanger designed to hang the heaviest apron on an apron rack or to hang it on a hook behind a door or wall hook. Raised end design prevents apron from slipping off the hanger.

Code	Details	Dimensions	Weight
ARHO	Apron Hanger (Open Loop)	19" W x 8"H	1lb
ARHC	Apron Hanger (Closed Loop)	19" W x 8"H	1lb



Vest & Skirt Apron Hanger

A chrome plated steel hanger with skirt hooks.

Code	Details	Dimensions	Weight
ARHO-SH	Vest & Skirt Hanger (Open Loop)	19" W x 8"H	1.6lb
ARHC-SH	Vest & Skirt Hanger (Closed Loop)	19" W x 8"H	1.6lb



Mobile Spider Apron Rack

An apron storage rack with eight spiral arms atop a post that is mounted on a star-shaped five leg base with nylon casters that is suitable for an Operating Room or a Cath Lab. Arms are slanted for better hold of aprons. This apron rack is available with optional glove holders, holding 3 pairs of gloves. Weight capacity 200lbs.

Code	Details	Dimensions	Weight
AR8	Mobile Spider Apron Rack	51" H x 24" W	30lbs
AR8-GH	Mobile Spider Apron Rack with Glove Holders	61 ^{1/2} " H x 24" W	36lbs

Mobile Spider Apron Rack




Top View of Rack with Glove Holders

Wall-Mounted Rack

for Half Aprons

Each rack accommodates 4 half-aprons. Half-aprons are not included. HAR-BV holds VELCRO® or Buckle half-aprons, and HAR-H holds Hoop half-aprons.



HAR-BV

HAR-H



Code	Details	Dimensions	Weight
HAR-BV	Wall-Mounted Rack only (Holds 4 VELCRO® or Buckle Half-Aprons)	21" L x 3½" H x 4" D	5.8/bs
HAR-H	Wall-Mounted Rack only (Holds 4 Hoop Half-Aprons)	17½" L x 4" H x 9¼" D	10.8lbs

Table Scatter Shield

We also manufacture custom-made screens and shields in all sizes and types. Supplied with or without hardware. See pages 22-24 for fabric options.



Code	Description
TDS-4824	Single Table Shield 48" x 24" - 0.5mm L.E.
TDS-4828	Single Table Shield 48" x 28" - 0.5mm L.E.
TDS-4836	Single Table Shield 48" x 36" - 0.5mm L.E.
TDS-3624	Single Table Shield 36" x 24" - 0.5mm L.E.
TDS-3628	Single Table Shield 36" x 28" - 0.5mm L.E.
TDS-3636	Single Table Shield 36" x 36" - 0.5mm L.E.
TDD-4824	Double Table Shield 48" x 24" - 0.5mm L.E.
TDD-4828	Double Table Shield 48" x 28" - 0.5mm L.E.
TDD-4836	Double Table Shield 48" x 36" - 0.5mm L.E.
TDD-3624	Double Table Shield 36" x 24" - 0.5mm L.E.
TDD-3628	Double Table Shield 36" x 28" - 0.5mm L.E.
TDD-3636	Double Table Shield 36" x 36" - 0.5mm L.E.
TD-CUST	Custom Table Shield
ATS-3228	Table Shield 32" x 28" Adjustable Height/Swing

Mobile Shield on T-Base

Adjustable height mobile shield that adjusts from 8" to 60". Excellent for patient and staff protection during imaging procedures. All steel T-Base provides excellent stability and confidence in no-tipping. Easy movement around the room and the unique base allows easy positioning around equipment and wall stands. Posi-Lock height adjustable arm with spring-loaded pin. Includes a durable nylon covered shield.

At lowest adjustment, will slide under a bench or table



Code	Description	Weight (approx)
MST	Frame Only	
MST-3024	Mobile Shield with Curtain 30" L x 24" W - 0.5mm L.E. with velcro attachment	27lbs
MST-3024-1	Mobile Shield with Curtain 30" L x 24" W - 1.0 mm L.E. with bolt attachment	37lbs

Mobile Shield on H-Base

Adjustable height: up to 60". Post Height: 60". Base Overall Dimensions: 24" W x 24" D. All steel H-base for maximum stability. Hardened chrome steel post. Posi-Lock height adjustable arm with spring-loaded pin. Includes a durable nylon covered shield with **Bolt attachment only**.

Code	Description	Weight (approx)
MSH	Frame Only	
MSH-3636	Mobile Shield with Curtain 36" x 36" - 0.5mm L.E. with bolt attachment	44lbs
MSH-4848	Mobile Shield with Curtain 48" x 48" - 0.5mm L.E. with bolt attachment	52lbs



Porta Shield

Adjustable mobile shield with soft curtain. Adjustable height from 24" to 48". Preset holes in the frame allow for swift and easy adjustment. Frame sits on four heavy duty casters for easy movement. Base measurement is 24" x 22".

Code	Description	Weight (approx)
PS5	Porta Shield with Curtain 24" x 36" - 0.5mm L.E.	40lbs
PS7	Porta Shield with Curtain 24" x 36" - 0.75mm L.E	42lbs
PS10	Porta Shield with Curtain 24" x 36" - 1.0mm L.E	45lbs
PS-CUST	Porta Shield - Custom	



Mobile Curtain Barrier

Excellent for patient protection during standing diagnostic procedures such as erect chest and cervical spine. Includes a durable nylon covered .5mm L.E. shield. Five castor star leg base. The base offers a sturdy 24-inch diameter, to ensure smooth and stable transportation around the room. It comes in three parts for quick and easy assembly. The adjustable telescoping hole allows the shield height to be set by the user. The top of the shield can be as low as 36 inches, or as high as 60 inches.

Code	Description	Weight (approx)
MS-2424	Mobile Shield with Curtain 24" x 24" - 0.5mm L.E. - with velcro attachment	23lbs
MS-2436	Mobile Shield with Curtain 24" x 36" - 0.5mm L.E. - with velcro attachment	25lbs
MS-2448	Mobile Shield with Curtain 24" x 48" - 0.5mm L.E. - with velcro attachment	28lbs
MS-CUST	Mobile Shield - Custom - 0.5mm L.E.	
MS2428SP	Mobile Shield with Split Panel Curtain 24" x 28" - 0.5mm L.E. - with bolt attachment	32lbs



Curtains

Curtains are ideal for areas where secondary or low-level radiation is present (Rooms and Doorways). Supplied with or without hard-ware. See pages 25-28 for fabric options.



Code	Description
DC5-27	Doorway Curtain 24"x 84" - 0.5mm L.E.
DC5-47	Doorway Curtain 48"x 84" - 0.5mm L.E.
DC10-27	Doorway Curtain 24"x 84" - 1.0mm L.E.
DC10-47	Doorway Curtain 48"x 84" - 1.0mm L.E.
DC-CUST	Doorway Curtain - Custom
RDC2	Room Divider Curtain / per sq. ft. - 0.25mm L.E.
RDC5	Room Divider Curtain / per sq. ft. - 0.5mm L.E.
RDC7	Room Divider Curtain / per sq. ft. - 0.75mm L.E.
RDC10	Room Divider Curtain / per sq. ft. - 1.0mm L.E.
TRA-4	Doorway / Room Divider 4' 0" Track
TRA-8	Doorway / Room Divider 8' 0" Track
TRO	Doorway / Room Divider Trolley (eacy)

Lead Vinyl Sheets

LeadX[®] is a smooth and flexible non-absorbent vinyl material that is easily cleaned.

Width options: 24" and 48".



0.25mm L.E.		0.5mm L.E.		1.0mm L.E.	
Code	Size	Code	Size	Code	Size
LX2-2424	24" x 24"	LX5-2424	24" x 24"	LX10-2424	24" x 24"
LX2-2436	24" x 36"	LX5-2436	24" x 36"	LX10-2436	24" x 36"
LX2-2448	24" x 48"	LX5-2448	24" x 48"	LX10-2448	24" x 48"
LX2-3648	36" x 48"	LX5-3648	36" x 48"	LX10-3648	36" x 48"
LX2-4848	48" x 48"	LX5-4848	48" x 48"	LX10-4848	48" x 48"
LX2-CUST	Custom	LX5-CUST	Custom	LX10-CUST	Custom
LX-PKG	Packaging Charge per Piece				