

UltraLite



- Stylish and comfortable wraparound frame.
- High-impact resistant plastic for safety.
- Colors: Black, Bronze, Carbon Fibre, Pink, Purple/Silver, Red, Silver and Yellow/Blue.
- Rx: Single and Bifocal.

Maxi



- Unisex shape suits any face and gender.
- Padded temple bars and nose bridge for extra comfort.
- The most durable product in our line.
- Colors: Black, Blue, Red and Silver.
- Rx: Single, Bifocal and Progressive.

Sport Wrap



- Solid, sturdy and extremely comfortable wraparound
- Built-in side-shields for safety.
- Colors: Black/Yellow and Brown/Orange.
- Rx: Single, Bifocal and Progressive.

Contour



- Designed to fit over existing prescription glasses.
- Black rubberized nylon frame.
- Lightweight and comfortable.
- Colors: Black and Tortoise.

Fitover - Regular



- Designed to fit over existing prescription glasses.
- Black rubberized nylon frame.
- Lightweight and comfortable.
- Colors: Black and Tortoise.

Fitover - Large



- Designed to fit over existing prescription glasses.
- Black rubberized nylon frame.
- Lightweight and comfortable.
- Colors: Black and Tortoise.

Racer



- Stylish wrap-around frame constructed of high impact TR-90 nylon.
- Rubberized soft grip temple bars.
- 8 base wrap around lenses for full protection around the orbit of the eye

The Lex



- Comfortable, lightweight and sturdy wide fit frame.
- Rubberized nose pads and temple bars provide a comfortable grip.
- Wraparound design ensures that your eyes are protected.



Webby Designer

- Very chic head-hugging design for ladies
- Durable nylon plastic high-impact safety frame.
- Chemical resistant.
- Colors: Black/White and Blue/Black.
- Single, Bifocal and Progressive.



Basic 52 Wrap

- Excellent peripheral vision with leaded glass side shields.
- Economical and lightweight with a spring hinge.
- Silicone, adjustable non-slip uni-fit bridge
- Rx: Single, Bifocal and Progressive.



Stars

- Premium protective women's eyewear.
- “Jackie-O” sunglasses style with large lenses.
- Excellent splash and side protection.
- Colors: Black/White and Red/Black.
- Rx: Single, Bifocal and Progressive.



Leopard

- Designer flair for style and performance.
- Head-hugging temple design with splash and side protection.
- Colors: Black and Tortoise.
- Rx: Single, Bifocal and Progressive.

EGM Mini Wraparound



- Stylish wrap-around frame constructed of high impact TR-90 nylon.
- Rubberized soft grip temple bars.
- Rubberized saddle bridge nose pad for extra comfort
- Stylish two-tone rubberized temple bars.

Splash Goggles - Regular



- Excellent peripheral vision with leaded glass side-shields.
- Economical and lightweight with a spring hinge.
- Comes with a silicone non-slip.
- Colors: As shown.
- Rx: Single, Bifocal and Progressive.

Splash Goggles - Large



- Excellent peripheral vision with leaded glass side-shields.
- Economical and lightweight with a spring hinge.
- Comes with a silicone non-slip.
- Colors: As shown.
- Rx: Single, Bifocal and Progressive.

Please Note: All glasses are 0.75mm L.E. Side-Shields (If equipped) are 0.50mm L.E. with the exception of models with integral .75mm lateral shield. All glasses include a protective case and strap.

Scatter Protector

(Full or Half Mask)

- Transparent ultralight lead acrylic.
- Shield is suspended from an adjustable headpiece.
- Both feature a headband support with an adjustable dial for circumference and angular positioning for that perfect fit and maximum protection.
- Non-Rx only.

CAUTION: Lead Equivalency = 0.10mm

Do not sterilize or use abrasive cleaning materials.

Full Mask: 8" Wrap-around design (Weight: 395g)

Half Mask: 5" Wrap-around design (Weight: 350g)

