



Open Hand Shield

- This shield resembles a paddle with a cuff. It is used to protect the top of the hand. One-size fits all with a lead-equivalency of 0.5mm.

Code	Description
GL-OH	Open Hand Shield



Open Palm Mitts

- These mitts are primarily designed for veterinary applications where fingers and thumbs are available to be used in a procedure. One-size fits all with a lead-equivalency of 0.5mm.

Code	Description
GL-OP	Open Palm Mitts



Open Palm with Thumb

- Identical to the Open-Palm Mitt with a closed thumb. One-size fits all with a lead-equivalency of 0.5mm.


Code	Description
GL-OPT	Open Palm Mitts with Thumb



Ducky Mitts

- These mitts are primarily designed for veterinary applications where fingers and thumbs must be protected. One-size fits all with a lead-equivalency of 0.5mm.

Code	Description
GL-DM	Ducky Mitts



Closed Mitts with Thumb

- Similar to the Ducky Mitt without the finger pouch. One-size fits all with a lead-equivalency of 0.5mm.

Code	Description
GL-DMT	Ducky Mitts with Thumb

Standard Colour: Royal Blue
Optional colors are available.

* Different colors shown are for deliniation purposes only.

Oven Style Mitts

An oven style mitt that offers 360° protection. One-size fits all with a lead-equivalency of 0.5mm.

Code	Description
GL-OM	Oven Style Mitts



Slit Mitts

A variation of the Oven Style Mitts with the upper section and thumb area open. One-size fits all with a lead-equivalency of 0.5mm.

Code	Description
GL-SM	Slit Mitts



15" Lead Vinyl Gloves

Each glove features an outer one-piece soft urethane molded design with an inner liner made from breathable nylon for comfort and is easily removed for inspection. Available individually or in pairs with a lead-equivalency of 0.5mm.

Code	Description
GL15	Seamless lead vinyl Glove - Pair
GL15L	Seamless lead vinyl Glove - Left
GL15R	Seamless lead vinyl Glove - Right



Molded One Piece Gloves

without an Inner Liner

These economical seamless lead vinyl gloves have no liner for easy washing inside and out. Size: Small, Medium and Large with a lead-equivalency of 0.5mm. LATEX FREE.

Code	Description
GL16-S	Molded One Piece Gloves - Small Pair
GL16-M	Molded One Piece Gloves - Medium Pair
GL16-L	Molded One Piece Gloves - Large Pair

